MANUAL DE INSTRUCCIONES / USER GUIDE
TABLERO ESPINAL / SPINE BOARD

REV. 2023/05

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01	MODELS.....	29
01.1	Products.....	29
01.2	Accessories.....	29
02	INTRODUCCIÓN.....	30
02.1	Using this guide.....	30
02.2	Legend of symbols.....	30
02.3	Service request.....	31
02.4	Life span.....	31
02.5	Demolition.....	31
02.6	Reference standards.....	31
02.7	Labelling.....	32
03	INTENDED USE.....	33
03.1	Intended use.....	33
03.2	Contraindications and side effects.....	33
03.3	Physical requirements of the operators.....	33
04	WARNINGS.....	34
04.1	General warnings.....	34
04.2	Specific warnings.....	36
04.3	Residual risk.....	37
05	DESCRIPTION OF THE PRODUCT.....	38
05.1	Main components.....	38
05.2	Technical data sheet.....	41
05.3	Main components (accessories).....	42
05.4	Technical data sheet (accessories).....	43
06	OPERATION.....	44
06.1	Transport and storage.....	44
06.2	Preparation.....	44
06.3	Functioning.....	45
07	GENERAL MAINTENANCE.....	47
07.1	Cleaning.....	47
07.2	Lubrication.....	47
07.3	Wear areas.....	47
07.4	Maintenance program.....	47
08	TRAINING REGISTER.....	49
09	MAINTENANCE REGISTER.....	50
10	LEGAL NOTICES.....	51
11	PRODUCT WARRANTY.....	52

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REV. 2023/05

01.1. Products

REFERENCE	UDI CODE	TRADE NAME
PA210-O	8437021314557	SPINE BOARD PA210 ORANGE
PA210-Y	8437021314564	SPINE BOARD PA210 YELLOW
PA212	8437021314571	PEDIATRIC SPINE BOARD
PA-252	8437021314588	ADULT AND PEDIATRIC SPINE BOARD
PA2400-00000	8437021314595	X-RAY SPINE BOARD
PA-06	8437021314601	X-RAY SPINE BOARD PA-06
PA-07	8437021314618	PEDIATRIC SPINE BOARD X-RAY
PA205	8436601870483	SHORT WOOD BOARD

01.2. Accessories

REFERENCE	UDI CODE	TRADE NAME
PA-130	8437021314526	SET OF 3 BELTS
80011	8437021314625	PEDIATRIC SPIDER STRAP FOR SPINE BOARD
PA-170	8437021314632	SPIDER STRAP FOR SPINE BOARD
PA-55	8437021314649	HEAD INMOBILIZER
PA-11	8437021314656	UNIVERSAL HEAD INMOBILISER
PA-241	8437021314663	PEDIATRIC HEAD INMOBILISER
PA-13	8436601870025	CPR BOARD PA-13

02.1. Using this guide



The manual provides instructions for use and maintenance of the product, as well as technical aspects, operation, spare parts and safety.

It is recommended before the operation of the product to read carefully this manual in order to avoid damages caused by a misuse, regardless of the level of experience acquired in the past with similar devices.

Do not lose this document. It should be accessible to any doubt that could appear by medical personnel. In case the device has some instructions for use of a different product than the one received, it is necessary to contact the manufacturer immediately before using the product.

Remember that a good use and maintenance are necessary for the proper operation of the product.

02.2. Legend of symbols

SYMBOL	EXPLANATION / DESCRIPTION
	MANUFACTURER symbol. This symbol is accompanied by the name and address of the manufacturer, adjacent to the symbol (PRODUCTOS METÁLICOS DEL BAGÉS S.L., Ctra. C-16 Km 59.5, 08650 Sallent (Barcelona)).
	Indicates the manufacturer's reference number to identify the medical device. PROMEBA, S.L. uses this symbol to set each internal reference for each configuration and business variant
	Indicates the manufacturer's serial number to identify a specific medical device.
	Indicates the manufacturing date. The symbol must be accompanied by a manufacturing date (yyyy-mm), adjacent to the symbol.
	It is placed to inform that the product is a "Medical Device"
	CE symbol without the intervention of a Notified Organism, as a Medical Device classified as Class I according to the EU Regulation 2017/745 on Medical Devices.
	Symbol for the Unique Device Identifier.
	Symbol "See instructions for use or operating instructions"
	Symbol "caution". This symbol is placed to warn of the need for the user to refer to important precautionary information in the operating instructions, such as warnings and cautions not otherwise found on the label.
	Symbol "Caution". For a general warning.
	Warning, crushing of hands
	Indicates the temperature limits. The upper and lower temperature limits should be indicated adjacent to the horizontal lines.

02.3. Service request

For information of the correct interpretation of the instruction manual, the use, maintenance, installation and restoration of the product, please contact Promeba customer service: T. 93 837 12 00, email: promeba@promeba.com or write to PROMEBA, S.L. - Ctra C-16 Km 59.5 · 08650 Sallent (Barcelona) · SPAIN.

Whenever you want to communicate any matter related to a product, you will be asked for the product reference (**REF**) and the serial number (**SN**). It is strongly recommended to keep the label containing this information and found on the original packaging of the product.

02.4. Life span

If used as described in the following instructions, this device has a useful life of 10 years from the date of purchase.

This useful life can be extended with annual reviews carried out by the manufacturer, which uses specialized and authorized internal and external technicians.

In case these annual checks are not carried out, the device must be disposed of according to the information in paragraph 02.5 and the manufacturer must be notified.

Only the manufacturer or an authorized center can extend the life of the device, if it meets the safety and legal requirements.

Promeba, S.L. will not accept any responsibility for malfunction or damage caused by the use of devices that have not been checked by the manufacturer or authorized center, or that have exceeded the maximum allowed useful life.

02.5. Demolition

When the devices are no longer suitable for use, if they have not been contaminated by any particular agent, they can be disposed of as normal solid waste, otherwise, follow the current demolition regulations.

02.6. Reference standards

REFERENCE	TITLE OF DOCUMENT
EU Regulation 2017/745	EU Regulation on Medical Devices

As a distributor or end user of the products manufactured and/or distributed by Promeba, S.L., you are strictly obliged to know the legal provisions in force in the country of destination of the goods, applicable to the devices to be supplied (including the regulations relating to technical specifications and/or safety requirements) and, therefore, understand the necessary requirements to ensure compliance of the products themselves with all the legal requirements of the territory.

02.7. Labelling

Each product incorporates an identification label, placed on the device itself and/or on the original packaging of the product. It must never be removed or covered. This label includes the data of the company which produced the device (🏭), the product code / reference (REF) of the product, the manufacturing date (📅), the serial number (SN), the UDI-PI code (UDI) and the CE distinctive (CE). Please keep these numbers so you can inform the dealer if necessary.

The UDI-PI code is made up of the following fields:

- (01)08437021314557: UDI code of the product.
- (11)230427: it gives information about the batch release date of the product following the format *aammdd*, where *aa* is the batch release year, *mm* is the month of the batch release and *dd* is the day of the batch release.
- (21)03241723033: identifies the serial number of the product giving information about the manufacturing order including the week when the product was produced.

03 INTENDED USE

REV. 2023/05

03.1. Intended use

The product SPINE BOARD is a device for lift and immobilize patients with possible spinal cord injuries and carry them to transport devices such as transport stretchers. With the use of the device, the patient can be picked up without changing his position and his trunk and extremities can be held up while being carried.

03.2. Contraindications and side effects

The contraindications resulted from the clinical analysis of the SPINE BOARD family of products are the following:

- Pain
- Compromised breathing caused by immobilization
- Pressure sores
- Need for additional radiological tests

03.3. Physical requirements of the operators

Promeba spine boards are destined to professional use only. It is expected that the devices will be used by first responders who have knowledge and experience with the immobilization and transport of patients with spinal cord injuries or various fractures.

The operators must be trained in efficient, effective and safe patient transport and must have the following minimum requirements:

- Physical capacity for operating the device.
- Be able to seize the device firmly with both hands.
- Have strong back, arms and legs for lifting, pushing and pulling the board.
- Have a good muscular coordination

It is recommended the employment of operators with strength, balance, coordination and common sense.

Patient loading procedures for extremely heavy patients, operations in rough terrain and in particular situations may require more operators.

The capacities of the various operators must be considered before determining his role in the employment of the board.

04.1. General warnings

1. The product must be used by trained personnel only, having attended specific training for this device and not for similar products.
2. Establish a maintenance and periodic control program and identify a reference manager. The person in charge of the ordinary maintenance of the device must guarantee the basic requirements provided by the manufacturer in these instructions for use. It is recommended that the subject check at least every 6 months if there are updated user manuals available or if there are any changes related to his / her product. This information is freely available on the website <http://promeba.com/>.
3. Training routines must be registered on a special register in which the names of those trained, of the trainers, date and place are indicated. This register, which will certify the eligibility of the operators to use the Promeba, S.L. device, has to be kept for a period of 10 years after the disposal of the device itself. This register will be made available to the competent authorities and/or manufacturer if requested.
4. Promeba, S.L. is always at your disposal to plan trainings on products.
5. Before carrying out any kind of operation on the appliance (training, installation, use), the operator must carefully read the enclosed instructions, paying particular attention to the correct safety precautions and to the procedures to be followed for installation and for correct use.
6. If the instructions belong to another device and not the device received, inform the manufacturer immediately and avoid use of the device.
7. In the case of any doubts as to the correct interpretation of the instructions, please contact Promeba, S.L. for any necessary clarifications.
8. Do not allow untrained persons to help during the use of the device, because they could cause damage to the patient or to themselves.
9. Before each use of device the perfect operating state of the device must be checked as specified in the instruction manual. If any damage or abnormalities which could in any way influence the correct functioning and the safety of the device, of the patient and/or of the user are detected, the device must be immediately removed from service and the manufacturer must be contacted.
10. If any failure or incorrect functioning of the device is detected, it must be immediately substituted with a similar item so that the rescue procedures are guaranteed without any interruption.
11. Use of the device in anyway other than described in this manual is forbidden.
12. Do not alter or modify in any way the device; any such interference could cause malfunctions and injury to the patient and/or rescuer.
13. The device must not in any way be tampered with (modification, adjustment, addition, replacement). In such cases all responsibility will be denied for any malfunctions or injuries caused by the device itself; moreover CE certification and product warranty will be considered void.

14. Ensure that all the necessary precautions are taken in order to avoid the hazards that can arise as the result of contact with blood or body fluids.
15. Register and store with these instructions: lot number, place and date of purchase, first date of use, date of checks, name of users and any other comments.
16. When the device is being used, the assistance of qualified staff must be guaranteed.
17. The product must not be exposed or come into contact with thermal sources of combustion and flammable agents, but must be stored in a dry, cool place, protected from light and sunlight.
18. Store and transport device in its original packaging.
19. Do not store the product under other more or less heavy materials, which can damage the structure of the product.
20. Position and adjust the device taking care not to cause any obstruction to rescuers and or any other rescue equipment.
21. Attention: laboratory testing, post production tests and instruction manuals cannot always consider every possible scenario for use. This means that in some cases the performance of the product could be notable different from results to date obtained. Instructions are continually being updated and are under tight surveillance of fully qualified staffs with adequate technical formation.
22. As a distributor or end users of products manufactured and/or marketed by Promebeba, S.L., you are strictly required to have a basic knowledge of any legal requirements existing in the goods final destination country applying to the devices contained in this supply (including laws and norms regarding technical specifications and/or safety requirements) and therefore you are also strictly required to have the necessary knowledge to guarantee all aspects regarding the total conformity of the products to the regulations in the relevant territory.
23. Promptly notify PROMEBEBA regarding any revisions to be made by manufacturer in order to guarantee the conformity of the product to the territory's legal specifications (including those resulting from rules and/or norms of other nature).
24. Act with all due care and diligence, and contribute to ensure conformity to general safety requirements of all devices marketed in the territory, by providing final users with all necessary information for carrying out periodical checks on their devices, as specified in the relevant user manual.
25. Actively contribute to product safety checks on products sold, by communicating any relevant risk analysis information both to the manufacturer and to any competent authorities so that the necessary actions can be promptly taken.
26. Be aware that in the event of any failure to conform to the above mentioned requirements you will be deemed fully responsible for all damages that might occur. Therefore we expressly disclaim any responsibility and/or liability for your non-compliance with the present "Regulatory provisions".

04.2. Specific warnings

1. Use only accessories/spare parts that are original or approved by Promeba, S.L. when performing any operation so as not to cause any alteration or modification to the device, otherwise we do not assume any responsibility for the correct operation or resulting damage of the device to the patient or the operator and the guarantee will be considered void, as well as compliance with the Regulation 2017/745.
2. Always respect the maximum capacity of the device, as indicated in this user's manual. Maximum load capacity means the total weight distributed according to the human anatomy. In determining the load of the total weight on the product, the operator must consider the weight of the patient, the equipment and the accessories. Moreover, the operator must consider that the overall dimensions of the patient do not reduce the functionality of the device.
3. Never leave the patient unassisted on the device, because he may be injured.
4. Do not use bleach to disinfect the device. Use a hydroalcoholic based disinfectant and wash with water.
5. The device and all its components, after washing, should be allowed to dry completely before storing.
6. Avoid contact with sharp objects.
7. Do not use the device if it is pierced, torn or frayed.
8. Make sure, before lifting, that the operators have a firm grip on the device.
9. Avoid pulling the device on rough surfaces.
10. The application of the device is not expected to last longer than the time required for first response operations and the subsequent phases of transportation to the nearest relief point. Do not use as a parking device.
11. To use the device, a minimum of two operators are needed in adequate physical conditions, so they must have strength, balance, coordination and common sense and must be trained in the device correct operation.
12. For particularly heavy patient loading and for rescue operations on steep terrain or in unusual circumstances, the presence of more operators is recommended.
13. Make sure the belts are properly fastened to the frame of the stretcher.
14. Do not operate in case the weight has not been distributed correctly.
15. Pay close attention to possible obstacles (water, ice, debris, etc.) in the path of the stretcher, as they could cause the operator to lose balance and compromise the proper functioning of the device. If you cannot clear the path of obstacles, choose an alternative path.
16. Condensation, water, ice and dust accumulations can affect the correct operation of the device, making it unpredictable and causing a sudden change in the weight that operators must carry.
17. Always immobilize the patient at least with the belts supplied by the manufacturer.

04.3. Residual risk

The residual risks listed below have been identified only with reference to the intended use of the device:

1. Use by untrained personnel may result in injury to patient, operator, or third parties.
2. Inappropriate disinfection procedures can create a risk of cross infection.
3. Failure to comply with the warnings for operators can create risks.
4. Failure to read and understand the product instructions can result in injury to the patient and operators.

05 PRODUCT DESCRIPTION *REV. 2023/05*

06.1. Main components

PA210-O / PA210-Y SPINE BOARD PA210 ORANGE / YELLOW

N°	DESCRIPTION OF THE COMPONENTS
1	Main spine board body
2	Handles and pins


PA212 PEDIATRIC SPINE BOARD

N°	DESCRIPTION OF THE COMPONENTS
1	Main spine board body
2	Handles and pins


PA-252 ADULT AND PEDIATRIC SPINE BOARD

N°	DESCRIPTION OF THE COMPONENTS
1	Main body of the adult spinal board
2	Main body of the pediatric spinal board
3	Pediatric spinal board fixation stops
4	Handles and elements for holding and fixing the adult spinal board
5	Handles and elements for holding and fixing the pediatric spinal board


PA2400-00000 X-RAY SPINE BOARD

N°	DESCRIPTION OF THE COMPONENTS
1	Main spine board body
2	Handles and pins


PA-06 X-RAY SPINE BOARD PA-06

N°	DESCRIPTION OF THE COMPONENTS
1	Main spine board body
2	Handles and pins


PA-07 PEDIATRIC SPINE BOARD X-RAY

N°	DESCRIPTION OF THE COMPONENTS
1	Main spine board body
2	Handles and pins

PA205 SHORT WOOD BOARD

N°	DESCRIPTION OF THE COMPONENTS
1	Main spine board body
2	Handles and pins
3	Immobilization belts

06.2. Technical data sheet


Models	PA210-O	PA210-Y	PA212	PA-252	PA-240	PA-06	PA-07	PA205
Length	1850mm	1850mm	1250mm	1830mm	1850mm	1850mm	1380mm	850mm
Width	425mm	425mm	425mm	450mm	340/410mm	410mm	350mm	400mm
Height	40mm	40mm	40mm	55mm	45mm	50mm	50mm	15mm
Weight	5,3kg	5,3kg	3,3kg	8kg	6,5kg	6,5kg	4kg	3,5kg
Max. Load	250kg	250kg	160kg	180kg	150kg	180kg	80kg	150kg
Color	Orange	Yellow	Orange	Black-Orange	Yellow	Yellow	Yellow	Brown
Material	ABS	ABS	ABS	HDPE	PE	PE	PE	Wood
Adult	X	X			X	X		X
Pediatric			X	X			X	

05 PRODUCT DESCRIPTION *REV. 2023/05*

06.3. Main components (accessories)


80011 PEDIATRIC SPIDER STRAP FOR SPINE BOARD
PA-170 SPIDER STRAP FOR SPINE BOARD

Nº	DESCRIPTION OF THE COMPONENTS
1	Cross tape with strap
2	Longitudinal tape
3	"V" tape

PA-55 HEAD IMMOBILIZER
PA-11 UNIVERSAL HEAD IMMOBILIZER
PA-241 PEDIATRIC HEAD IMMOBILIZER

Nº	DESCRIPTION OF THE COMPONENTS
1	Cushion
2	Base
3	Belts

PA-13 CPR BOARD PA-13

N°	DESCRIPTION OF THE COMPONENTS
1	Main spine board body
2	Handles and pins

06.4. Technical data sheet (accessories)

Models	80011	PA-170	PA-55	PA-11	PA-241	PA-13
Length			420mm	420mm	260mm	640mm
Width			260mm	260mm	340mm	430mm
Height			180mm	180mm	140mm	80mm
Weight	0,5kg	0,8kg	0,8kg	0,8kg	0,5kg	1,5kg
Max. Load	250kg	250kg				
Color	Orange			Red	Blue	Blue
Material	Polypropylene / Nylon	Polypropylene / Nylon	PU Foam + PVC	PU Foam + PVC	PU Foam + PVC	PE
Adult		X	X	X		X
Pediatric	X				X	

06.1. Transport and storage

Before transporting the device, make sure that it is correctly packaged.

During transport always fix the load. If piling up is necessary always follow the scheme shown on figure 1.



Transport the leveled load and following all precepts and rules for the transport of loads, ensuring also that there are no risks of shocks, bumps or falls during the transport itself.

To unpack the device place the box on a flat, stable surface and carefully open the seal. Remove the device from the inside of the box following the scheme shown on figure 2. Keep the original packaging for use in case of any further transport and for storage.

Damage to the device caused during transport and handling is not covered by the guarantee. Repairs or replacement of the damaged parts are the responsibility of the client.

The device must be stored in a dry, cool area away from direct sunlight. It must not be placed in contact with any substances or chemical agents which could cause damage and reduce safety characteristics.

The storage temperature must be between -20°C to 40°C.

The product must not be manipulated in any way (modification, adjustment, addition, substitution). In such cases, all responsibility for any malfunction or injury caused by the product itself will be disclaimed; furthermore, CE certification and product warranty will be considered void.

06.2. Preparation

On receipt of the product:

- Remove the packaging and display the material so that all components are visible.
- Check that all the components/pieces on the accompanying list are present.

The appliance must be checked before every use so as to reveal any working abnormalities and/or damage caused by transport and/or storage. In particular, check:

- General functionality of the device.
- Cleanliness of the device (remember that the failure of cleaning may cause the risk of cross infections).
- Absence of cuts, holes, tears on the structure, including the straps.
- Integrity of components.

If the above conditions are met, the device may be considered ready for use; otherwise it must be immediately removed from service and the manufacturer must be contacted.

06.3. Functioning


Guidelines for using the spine board

Before using the spine board, you must carefully read this manual and the manual for all accessories that are intended to be used with the board, such as belts, spider belts or head immobilizers.

You must strictly follow the guidelines of your emergency department and/or local and regional guidelines before placing the patient on the spinal board.

Carefully assess the need for use of a cervical collar, spider belt, head immobilizer, or other means of immobilization in addition to restraint belts.

The PA-252 tabletop is equipped with 4 zones of different depths in relation to the patient's support surface. The correct position in each service must be evaluated by the lifeguard / rescue professional, to improve the alignment of the cervical tract, compensating for the different occipital prominence typical of pediatric patients.

To separate the pediatric part of the adult spine board (in the PA-252 model), turn the red stops 90° to be able to remove the pediatric board.

When you have finished using it, place it back inside the adult board, turning the same stops on which you had previously turned, making sure that the two devices are correctly joined, as shown in the image.

Spinal boards can be used to make preliminary diagnoses to X-rays, to confirm or refute the hypothesis of the presence of spinal lesions in the patient. Still, these devices should not be used in MRI settings.

Immobilization of the patient:

After positioning the patient, it is essential to apply the support belts or spider belts according to the procedure established for these devices.

It will be the lifeguard / rescue professional who must decide if the patient's condition requires or allows it to apply the head, ankle or cervical collar immobilizer until ensuring adequate and safe immobilization for transport.

Once the patient is immobilized, the device can be lifted and transported. Never use cranes, slings or other means to lift the device.


Lifting is only allowed manually and by an adequate number of operators based on the patient's weight. (Consider also the weight of the device and any accessories that may be present) The experience gathered by Promeba advises to carry out the lifting with four (4) technicians if possible, although in some cases two (2) would be sufficient.

06 OPERATION

REV. 2023/05

Operator positioning should be symmetrical and allow safe and stable lifting of both the head and foot sides of the spinal board.

The image below shows a placement recommendation from EMTs to make the lift as safe and stable as possible.

When there are only two operators, it is recommended that they be located either in a central position or somewhat displaced towards the part of the patient's head, which is where the greatest weight is concentrated (round in green).

When the elevation can be done with four technicians, one must be placed at each end of the board as the red circles are located in the image above. Finally, for cases in which the weight of the patient recommends more operators, the recommended location is as the rounds are distributed around the board.

The maximum load for each operator must never exceed the maximum load allowed by occupational safety regulations.

07 GENERAL MAINTENANCE REV. 2023/05**07.1. Cleaning**

It is essential to keep the equipment clean to ensure proper use and durability of the set. Thorough cleaning should be carried out periodically, especially in areas exposed to dirt that can accumulate dirt such as corners or holes.

Do not use high pressure cleaning systems or bleach to disinfect the product, they can damage it. Instead, use a hydroalcoholic-based disinfectant, wash it with water and let it dry naturally, do not use direct heat sources to dry.

Non-compliance with cleaning operations can lead to the risk of cross infections due to the presence of secretions and/or residues. During all control and disinfection operations, the operator must wear appropriate personal protective equipment, such as gloves, glasses, etc.

When the cleaning procedure includes a disinfectant product, carefully follow the instructions of the manufacturer of the product used regarding the method of application and contact time. Make sure that you have taken all the necessary precautions to ensure that there are no risks of cross-infection or contamination of patients and operators.

07.2. Lubricating

Due to its intrinsic characteristics, the spinal board family of products does not require lubrication, since it does not have moving parts.

07.3. Wear areas

Periodic inspection of the different components of the system for signs of wear is a preventative measure that can reduce breakdowns. Check especially belts, straps, textiles and fastening handles as points with more wear.

In the event that a worn component needs to be changed, you should contact our sales department for more information on ordering spare parts and their installation.

The repair or replacement of product components must necessarily be carried out by qualified personnel, that is, who have received training from the Manufacturer, Promebe S.L. and original spare parts should always be used.

Promebe S.L. declines any responsibility for any damage, direct or indirect, resulting from improper use of spare parts and/or any repair intervention carried out by unauthorized persons.

07.4. Maintenance program

Although the frequency of the controls is determined by factors such as legal requirements, the type of use, the frequency of use, the environmental conditions during use and storage, it is placed a recommendation on the maintenance that the equipment must have for correct operation throughout its useful life.

All maintenance and revision activities must be recorded and documented with the corresponding technical intervention reports, and the documentation must be kept for at least

07 GENERAL MAINTENANCE *REV. 2023/05*

10 years from the end of the product's useful life and must be made available to the competent Authorities and/or from the Manufacturer, when requested.

Maintenance review	Every use	When necessary	Every month	Every year
Disinfect	X		X	X
Clean		X	X	X
Inspect	X	X	X	X

For both the disinfection and cleaning procedures, point 7.1 of this manual must be followed.

During the inspection procedure, which must be done after each use, the following tasks must be carried out:

- Check that all components are present.
- Check the integrity of the device, that is, that there are no breaks, cracks, holes or cuts.
- Check the wear of the product following section 7.3 of this user guide.
- Check that the moving parts slide correctly.
- Check that the accessories provided are present, intact and that they function correctly.

10 LEGAL NOTICES

REV. 2023/05

- This document may contain technical inaccuracies or typographical mistakes.
- Changes are periodically added to the information herein; these changes will be incorporated in new editions of the publication.
- Promeba, S.L. reserves the right to make any modification or improvement in the products described in this publication if it is appropriate.
- Promeba, S.L. may have patents or patent applications which address themes described in this document. The possession of this document does not entitle no license to these patents. The information contained in this document does not affect or change the specifications or product specifications or warranties of Promeba, S.L.
- No part of this document shall operate as express or implied license or compensation under the intellectual property rights of Promeba, S.L. or .third party.
- All the information contained in this document has been obtained in specific environments and it is presented as an illustration. The results obtained in other operating environments may vary.
- Promeba, S.L. can use or distribute the information that the customer provides in whatever way they see fit, without incurring any obligation with the client.

11 PRODUCT WARRANTY *REV. 2023/05*

Promeba, S.L. guarantees that its products have satisfactorily passed all established quality controls, both functional and material. The duration of the guarantee is 2 years from the date of purchase of the product.

This guarantee will be valid only when it is presented with the original invoice or proof of purchase (indicating the date of purchase, model and the name of the distributor) together with the defective product during the period covered by the guarantee. Promeba, S.L. reserves the right not to offer the free warranty service if the indicated documents are not presented or if the information they contain is incomplete or illegible.

1. This warranty will not apply if the model name or serial number of the product has been altered, erased, disappeared or becomes illegible.

2. This guarantee does not cover transportation costs or risks derived from the transportation of your product to and from Promeba, S.L.

3. This warranty does not cover any of the following:

- a) Periodic maintenance and repair or replacement of parts derived from normal wear and tear.
- b) Expendable material (components that are expected to need periodic replacement during the life of the product, such as belts, straps, etc.).
- c) Damages or defects derived from the improper use, operation or treatment of the product and not due to normal use of the product.
- d) Damages derived from:
 - i. Misuse, including:
 - Treatment that results in damages or physical, superficial or appearance changes of the product.
 - Installation, use or storage of the product in a way that does not

respect the instructions described by Promeba, S.L.

- Maintenance of the product in a way that does not respect the instructions of Promeba, S.L. for proper maintenance.

- Installation or use of the product in a way that does not respect the technical or safety regulations of the country where it is used or installed.

- ii. Use of components not provided with the product or incorrect installation of accessory parts not previously tested.

- iii. States or defects of the system in which the product is used or incorporated with the exception of other products Promeba, S.L. designed for use with the product.

- iv. Use of the product with accessories, peripheral units and other products of a type, condition or standards not established by Promeba, S.L.

- v. The manufacturer or distributor will be solely responsible for deciding whether to send the parts for repair, or to replace the product in its entirety. In no case will operators be sent to repair or replacement of the product.

Except in the cases mentioned above, Promeba, S.L. will make no warranties in relation to product, performance, accuracy, reliability, or fitness for logic purpose of the equipment or other type. If this exception is not legal or contemplated by current law, Promeba, S.L. will limit or exclude your warranties only to the extent permitted by applicable law.

The only obligation on the part of Promeba, S.L. in connection with this warranty is to repair or replace parts subject to the terms and conditions of this warranty.

Promeba, S.L. is not responsible for loss or damage to products, this warranty or others, including economic loss or non-assessable

11 PRODUCT WARRANTY *REV. 2023/05*

damage; the price paid for the product; loss of profits, income, information, usufruct or use of the product or associated products or indirect, accidental or critical loss or damage.

This clause refers to whether the loss or damage is due to deterioration or inoperability of the associated product due to defects or unavailability of Promeba, S.L., which has caused downtime, loss of user time or business interruption.

In cases where the law prohibits or limits these liability exclusions, Promeba, S.L. will exclude or limit his liability only to the extent permitted by applicable law.

For example, there are countries that prohibit the exclusion or limitation of damages caused by negligence, reckless negligence, willful misconduct, fraud and similar actions. The responsibility of Promeba, S.L. In this guarantee will not exceed, in any case, the price paid for the product, but if the current law allows only limitations of greater responsibilities, these will apply.

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*All rights reserved. Variations can be done whitout notice.
Promeba, S.L. is to be considered not responsible for damages caused by the lack or the wrongness of the information here mentioned.*

SCAN QR CODE

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Ctra C-16 Km 59.5 · 08650 Sallent (Barcelona) · SPAIN

T. 93 837 12 00

promeba@promeba.com · www.promeba.com