**MANUAL DE INSTRUCCIONES / USER GUIDE**  
**RAÍL SERIE PB-420 / RAIL SERIE PB-420**

REV. 2024/12

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01	MODELS.....	43
01.1	Products.....	43
01.2	Accessories.....	43
02	INTRODUCCIÓN.....	44
02.1	Using this guide.....	44
02.2	Legend of symbols.....	44
02.3	Service request.....	45
02.4	Life span.....	45
02.5	Demolition.....	45
02.6	Reference standards.....	45
02.7	Labelling.....	46
03	INTENDED USE.....	47
03.1	Intended use.....	47
03.2	Contraindications and side effects.....	47
03.3	Physical requirements of the operators.....	47
04	WARNINGS.....	48
04.1	General warnings.....	48
04.2	Specific warnings.....	50
04.3	Residual risk.....	51
05	DESCRIPTION OF THE PRODUCT.....	52
05.1	Main components.....	52
05.2	Technical data sheet.....	60
05.3	Main components (accessories).....	61
05.4	Technical data sheet (accessories).....	61
06	OPERATION.....	62
06.1	Transport and storage.....	62
06.2	Preparation.....	62
06.3	Functioning.....	64
07	GENERAL MAINTENANCE.....	67
07.1	Cleaning.....	67
07.2	Lubrication.....	67
07.3	Wear areas.....	67
07.4	Maintenance program.....	68
08	SPARE PARTS.....	69
09	TRAINING REGISTER.....	75
10	MAINTENANCE REGISTER.....	76
11	LEGAL NOTICES.....	77
12	PRODUCT WARRANTY.....	78

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REV. 2024/12

**01.1. Products**

REFERENCE	UDI CODE	TRADE NAME
PB4200-00000	8437001519781	STRETCHER SUPPORT WITH DISPLACEMENT MOVEMENT
PB4201-00000	8437001519804	LATERAL MOVEMENT BASIC STRETCHER SUPPORT PB-4201
PB4202-00000	8437001519828	MODULAR BASIC STRETCHER SUPPORT PB420-2
PB4203-00000	8437001519842	BASIC PLANE STRETCHER SUPPORT PB-420-3
PB4204-00000	8437001519859	DOUBLE BASIC STRETCHER SUPPORT PB420-4
PB4206-00000	8437001519866	STRETCHER SUPPORT + LATERAL DISP. PB420-6
PB4207-00000	8437001519873	PB420-2 STRETCHER SUPPORT WITH LOCK
PB4230-00000	8437001519897	STRETCHER SUPPORT PB-423 AUTOMATIC OPEN
PB4231-00000	8436601870438	STRETCHER SUPPORT SET ALUMINUM RAL.5010
PB4210-00000	8436601870445	STRETCHER SUPPORT WITH DRAWER, ANCHORAGE AND RAMP

**01.2. Accessories**

REFERENCE	UDI CODE	TRADE NAME
PB4200-03172	8437001519996	SHORT RAMP SET PROMEBA
PB4200-03072	8437021314021	PROMEBA LONG RAMP SET – STD
PB4200-03310	8437021314069	LONG RAMP SET - 450MM
PB4200-03271	8437021314076	LONG RAMP SET - 500mm
PB4200-03290	8437021314090	LONG RAMP SET - 550mm
PB4205-03090	8437021314113	ULTRA LONG RAMP SET
PB4100-05150	8437021314120	LONG RAMP SET 900mm

## 02 INTRODUCTION

REV. 2024/12

### 02.1. Using this guide

The manual provides instructions for use and maintenance of the product, as well as technical aspects, operation, spare parts and safety.



It is recommended before the operation of the product to read carefully this manual in order to avoid damages caused by a misuse, regardless of the level of experience acquired in the past with similar devices.

Do not lose this document. It should be accessible to any doubt that could appear by medical personnel. In case the device has some instructions for use of a different product than the one received, it is necessary to contact the manufacturer immediately before using the product.

Remember that a good use and maintenance are necessary for the proper operation of the product.

Each product includes an identifying adhesive label with the serial number and model. Keep these numbers so that you can refer them to the dealer if necessary.

### 02.2. Legend of symbols

SYMBOL	EXPLANATION / DESCRIPTION
	MANUFACTURER symbol. This symbol is accompanied by the name and address of the manufacturer, adjacent to the symbol (PRODUCTOS METÁLICOS DEL BAGÉS S.L., Ctra. C-16 Km 59.5, 08650 Sallent (Barcelona)).
	Indicates the manufacturer's reference number to identify the medical device. PROMEBA, S.L. uses this symbol to set each internal reference for each configuration and business variant
	Indicates the manufacturer's serial number to identify a specific medical device.
	Indicates the manufacturing date. The symbol must be accompanied by a manufacturing date (yyyy-mm), adjacent to the symbol.
	It is placed to inform that the product is a "Medical Device"
	CE symbol without the intervention of a Notified Organism, as a Medical Device classified as Class I according to the EU Regulation 2017/745 on Medical Devices.
	Symbol for the Unique Device Identifier.
	Symbol "See instructions for use or operating instructions"
	Symbol "caution". This symbol is placed to warn of the need for the user to refer to important precautionary information in the operating instructions, such as warnings and cautions not otherwise found on the label.
	Symbol "Caution". For a general warning.
	Warning, crushing of hands
	Indicates the temperature limits. The upper and lower temperature limits should be indicated adjacent to the horizontal lines.

**02 INTRODUCTION**

REV. 2024/12

**02.3. Service request**

For information of the correct interpretation of the instruction manual, the use, maintenance, installation and restoration of the product, please contact Promeba customer service: T. 93 837 12 00, email: promeba@promeba.com or write to PROMEBA, S.L. - Ctra C-16 Km 59.5 - 08650 Sallent (Barcelona) - SPAIN.

Whenever you want to communicate any matter related to a product, you will be asked for the product reference (**REF**) and the serial number (**SN**). It is strongly recommended to keep the label containing this information and found on the original packaging of the product.

**02.4. Life span**

**If used as described in the following instructions, this device has a useful life of 10 years from the date of purchase.**

This useful life can be extended with annual reviews carried out by the manufacturer, which uses specialized and authorized internal and external technicians.

In case these annual checks are not carried out, the device must be disposed of according to the information in paragraph 02.5 and the manufacturer must be notified.

Only the manufacturer or an authorized center can extend the life of the device, if it meets the safety and legal requirements.

Promeba, S.L. will not accept any responsibility for malfunction or damage caused by the use of devices that have not been checked by the manufacturer or authorized center, or that have exceeded the maximum allowed useful life.

**02.5. Demolition**

When the devices are no longer suitable for use, if they have not been contaminated by any particular agent, they can be disposed of as normal solid waste, otherwise, follow the current demolition regulations.

**02.6. Reference standards**

REFERENCE	TITLE OF DOCUMENT
EU Regulation 2017/745	EU Regulation on Medical Devices
EN 1789:2007+A2:2014	Medical vehicles and their equipment - Road ambulances
UNE-EN 1865-5:2012	Patient handling equipment used in road ambulances - Part 5: Stretcher support


As a distributor or end user of the products manufactured and/or distributed by Promeba, S.L., you are strictly obliged to know the legal provisions in force in the country of destination of the goods, applicable to the devices to be supplied (including the regulations relating to technical specifications and/or safety requirements) and, therefore, understand the necessary requirements to ensure compliance of the products themselves with all the legal requirements of the territory.

**02 INTRODUCTION**

REV. 2024/12

**02.7. Labelling**

Each product incorporates an identification label, placed on the device itself and/or on the original packaging of the product. It must never be removed or covered. This label includes the data of the company which produced the device (🏭), the product code / reference (REF) of the product, the manufacturing date (📅), the serial number (SN), the UDI-PI code (UDI) and the CE distinctive (CE). Please keep these numbers so you can inform the dealer if necessary.

The UDI-PI code is made up of the following fields:

- (01)08437001519781: UDI code of the product.
- (11)230805: it gives information about the batch release date of the product following the format *aammdd*, where *aa* is the batch release year, *mm* is the month of the batch release and *dd* is the day of the batch release.
- (21)01233123001: identifies the serial number of the product giving information about the manufacturing order (0123) including the week when the product was produced (3123) and a correlative number depending on the units that has been produced in that production order (001).

## **03 INTENDED USE**

REV. 2024/12

### **03.1. Intended use**

The product BASIC STRETCHER SUPPORT is a device indicated to load and unload transportation stretchers from an emergency vehicle transport while supporting them, remaining anchored during the movement of the vehicle.

### **03.2. Contraindications and side effects**

Any contraindication or side effect have been resulted from the clinical analysis of the BASIC STRETCHER SUPPORT family of products.

### **03.3. Physical requirements of the operators**

Promeba rails are intended for professional use only. The device is intended to be used by medical transport technicians or ambulance technicians who have knowledge and experience with patient transport.

Operators must be able to assess the integrity of the fixing, knowing the mechanism for fixing the stretcher to the support and being able to ensure its correct operation. They must also be able to assess any anomaly in the system/assembly by communicating the problem to those responsible and stopping the use of the devices if necessary to ensure patient safety.

The capacities of the various operators must be considered before determining his role in the employment of the board.

**04.1. General warnings**

1. The product must be used by trained personnel only, having attended specific training for this device and not for similar products.
2. Establish a maintenance and periodic control program and identify a reference manager. The person in charge of the ordinary maintenance of the device must guarantee the basic requirements provided by the manufacturer in these instructions for use. It is recommended that the subject check at least every 6 months if there are updated user manuals available or if there are any changes related to his / her product. This information is freely available on the website <http://promeba.com/>.
3. Training routines must be registered on a special register in which the names of those trained, of the trainers, date and place are indicated. This register, which will certify the eligibility of the operators to use the Promeba, S.L. device, has to be kept for a period of 10 years after the disposal of the device itself. This register will be made available to the competent authorities and/or manufacturer if requested.
4. Promeba, S.L. is always at your disposal to plan trainings on products.
5. Before carrying out any kind of operation on the appliance (training, installation, use), the operator must carefully read the enclosed instructions, paying particular attention to the correct safety precautions and to the procedures to be followed for installation and for correct use.
6. If the instructions belong to another device and not the device received, inform the manufacturer immediately and avoid use of the device.
7. In the case of any doubts as to the correct interpretation of the instructions, please contact Promeba, S.L. for any necessary clarifications.
8. Do not allow untrained persons to help during the use of the device, because they could cause damage to the patient or to themselves.
9. Before each use of device the perfect operating state of the device must be checked as specified in the instruction manual. If any damage or abnormalities which could in any way influence the correct functioning and the safety of the device, of the patient and/or of the user are detected, the device must be immediately removed from service and the manufacturer must be contacted.
10. If any failure or incorrect functioning of the device is detected, it must be immediately substituted with a similar item so that the rescue procedures are guaranteed without any interruption.
11. Use of the device in anyway other than described in this manual is forbidden.
12. Do not alter or modify in any way the device; any such interference could cause malfunctions and injury to the patient and/or rescuer.
13. The device must not in any way be tampered with (modification, adjustment, addition, replacement). In such cases all responsibility will be denied for any malfunctions or injuries caused by the device itself; moreover CE certification and product warranty will be considered void.

## **04 WARNINGS**

*REV. 2024/12*

14. Use only original or PROMEBA approved components/spare parts and/or accessories to carry out all operations without causing alterations or modifications to the device; Otherwise, any responsibility for the operation is declined.
15. Ensure that all the necessary precautions are taken in order to avoid the hazards that can arise as the result of contact with blood or body fluids.
16. Register and store with these instructions: lot number, place and date of purchase, first date of use, date of checks, name of users and any other comments.
17. When the device is being used, the assistance of qualified staff must be guaranteed.
18. The product must not be exposed or come into contact with thermal sources of combustion and flammable agents, but must be stored in a dry, cool place, protected from light and sunlight.
19. Store and transport device in its original packaging.
20. Do not store the product under other more or less heavy materials, which can damage the structure of the product.
21. Attention: laboratory testing, post production tests and instruction manuals cannot always consider every possible scenario for use. This means that in some cases the performance of the product could be notable different from results to date obtained. Instructions are continually being updated and are under tight surveillance of fully qualified staffs with adequate technical formation.
22. As a distributor or end users of products manufactured and/or marketed by Promebeba, S.L., you are strictly required to have a basic knowledge of any legal requirements existing in the goods final destination country applying to the devices contained in this supply (including laws and norms regarding technical specifications and/or safety requirements) and therefore you are also strictly required to have the necessary knowledge to guarantee all aspects regarding the total conformity of the products to the regulations in the relevant territory.
23. Promptly notify PROMEBA regarding any revisions to be made by manufacturer in order to guarantee the conformity of the product to the territory's legal specifications (including those resulting from rules and/or norms of other nature).
24. Act with all due care and diligence, and contribute to ensure conformity to general safety requirements of all devices marketed in the territory, by providing final users with all necessary information for carrying out periodical checks on their devices, as specified in the relevant user manual.
25. Actively contribute to product safety checks on products sold, by communicating any relevant risk analysis information both to the manufacturer and to any competent authorities so that the necessary actions can be promptly taken.
26. Be aware that in the event of any failure to conform to the above mentioned requirements you will be deemed fully responsible for all damages that might occur. Therefore we expressly disclaim any responsibility and/or liability for your non-compliance with the present "Regulatory provisions".
27. Follow internal procedures and protocols approved by your organization.

28. Follow approved Emergency Medical Service procedures for patient immobilization and transportation.

### **04.2. Specific warnings**

1. Use only accessories/spare parts that are original or approved by Promeba, S.L. when performing any operation so as not to cause any alteration or modification to the device, otherwise we do not assume any responsibility for the correct operation or resulting damage of the device to the patient or the operator and the guarantee will be considered void, as well as compliance with the Regulation 2017/745.
2. Always respect the maximum capacity of the device, as indicated in this user's manual. Maximum load capacity means the total weight distributed according to the human anatomy. In determining the load of the total weight on the product, the operator must consider the weight of the patient, the equipment and the accessories. Moreover, the operator must consider that the overall dimensions of the patient do not reduce the functionality of the device.
3. Inside the ambulance, place and adjust the device taking care not to obstruct other parts or other equipment. Also consider lateral displacement if it is available in your unit.
4. Never leave the patient unassisted on the device, because he may be injured.
5. Do not use bleach to disinfect the device. Use a hydroalcoholic based disinfectant and wash with water.
6. The device and all its components, after washing, should be allowed to dry completely before storing.
7. Avoid contact with sharp objects.
8. Do not use the device if it is pierced, deformed or shows signs of loss of integrity.
9. Avoid dragging rough objects on the surface of the device.
10. Make sure the spine board or scoop stretcher is securely secured by using the rail belts before moving the ambulance. If it is not properly fixed, there is a risk of impact.
11. Do not operate in case the weight has not been distributed correctly.
12. Condensation, water, ice and dust accumulations can affect the correct operation of the device, making it unpredictable.
13. Always immobilize the patient at least with the belts supplied by the manufacturer during any movement of the stretcher, including loading and unloading to the ambulance.
14. Do not use the device as a support for climbing or leaning on and reaching high points in the ambulance. The product may be damaged and the technician or patient may be harmed.
15. Disinfection activities must be carried out in accordance with the validated cycle parameters, indicated in the specific technical standards.

## **04 WARNINGS**

REV. 2024/12

16. The stretcher used together with the stretcher support, must be installed in accordance with the specifications provided by the manufacturer and ensuring its compatibility with technical specifications of the stretcher support.

17. Promeba S.L. shall not be liable if unauthorized / incompatible stretcher fastener is used.

18. The installation of the stretcher support on the vehicle floor must be carried out ensuring compliance with the safety requirements of EN 1789. In order to be in compliance with the requirements of EN 1789, the installer has the duty to check the fixation based on the exact position on the floor and based on the structure of the floor itself.

19. Do not place body parts under the stretcher support; this may cause crushing risks when the stretcher support is moved. Never put your feet or hands under the device.

20. Installation carried out by untrained personnel, could result in inadequate fastening of the device inside the vehicle, with consequent risks related to its instability or mobility.

### **04.3. Residual risk**

The residual risks listed below have been identified only with reference to the intended use of the device:



1. Do not place body parts under the stretcher support; this may cause crushing risks when the stretcher support is moved. Never put your feet or hands under the device.

2. Installation carried out by untrained personnel, could result in inadequate fastening of the device inside the vehicle, with consequent risks related to its instability or mobility. 1. Use by untrained personnel may result in injury to patient, operator, or third parties.

**05.1. Main components**



**PB4200-00000 STRETCHER SUPPORT WITH  
DISPLACEMENT MOVEMENT**

<b>N°</b>	<b>DESCRIPTION OF THE COMPONENTS</b>
1	Stainless steel ramp
2	Accommodation for scoop stretcher
3	Ramp inclination regulator
4	Lateral displacement base
5	Lateral displacement pedal


**PB4201-00000 LATERAL MOVEMENT BASIC STRETCHER  
SUPPORT PB-4201**

<b>N°</b>	<b>DESCRIPTION OF THE COMPONENTS</b>
1	Stainless steel ramp
2	Lateral displacement base
3	Lateral displacement pedal


**PB4202-00000 MODULAR BASIC STRETCHER SUPPORT PB420-2**

<b>N°</b>	<b>DESCRIPTION OF THE COMPONENTS</b>
1	Stainless steel ramp
2	Accommodation for scoop stretcher
3	Ramp inclination regulator

**PB4203-00000 BASIC PLANE STRETCHER SUPPORT PB-420-3**

N°	DESCRIPTION OF THE COMPONENTS
1	Stainless steel ramp


PB4204-00000 DOUBLE BASIC STRETCHER SUPPORT PB420-4

N°	DESCRIPTION OF THE COMPONENTS
1	Stainless steel ramp
2	Accommodation for scoop stretcher
3	Accommodation for scoop stretcher
4	Ramp inclination regulator





PB4206-00000 STRETCHER SUPPORT + LATERAL DISP. PB420-6

N°	DESCRIPTION OF THE COMPONENTS
1	Stainless steel ramp
2	Accommodation for scoop stretcher
3	Accommodation for scoop stretcher
4	Ramp inclination regulator
5	Lateral displacement base
6	Lateral displacement pedal


**PB4210-00000 STRETCHER SUPPORT WITH DRAWER, ANCHORAGE AND RAMP**

N°	DESCRIPTION OF THE COMPONENTS
1	Stainless steel ramp
2	Accommodation for scoop stretcher
3	Lock for the stretcher included in the support

**PB4231-00000 STRETCHER SUPPORT SET ALUMINUM RAL.5010**

N°	DESCRIPTION OF THE COMPONENTS
1	Aluminium ramp
2	Automatic lock included in the rail


**05.2. Technical specifications / data sheet**

<b>Models</b>	<b>PB-420/0</b>	<b>PB420/1</b>	<b>PB-420/2</b>	<b>PB-420/3</b>	<b>PB-420/4</b>	<b>PB-420/6</b>	<b>PB-420/7</b>	<b>PB-423/0</b>	<b>PB-423/1</b>	<b>PB-421/0</b>
Length	2090mm	2060mm	2025mm	2010mm	2025mm	2080mm	2025mm	1910mm	1910mm	1904mm
Width	600mm	600mm	600mm	600mm	600mm	600mm	600mm	640mm	640mm	625mm
Height	170mm	70mm	100mm	55mm	200mm	270mm	100mm	125mm	125mm	129mm
Weight	45kg	40kg	41,5kg	30kg	60kg	65kg	41,5kg	30kg	30kg	30kg
Max. Load	350kg	350kg	380kg	380kg	380kg	350kg	380kg	250kg	380kg	380kg
Color	Blue	Blue	Blue	Blue	Blue	Blue	Blue I	Grey	Blue	Blue
Material	Steel	Steel	Steel	Steel	Steel	Steel	Steel	Alu	Alu	Steel
Lat. Disp.	300mm	300mm	NO	NO	NO	300mm	NO	NO	NO	NO
Lock	Not inc.	Not inc.	Not inc.	Not inc.	Not inc.	Not inc.	Inc.	Inc.	Inc.	Inc.

**05.3. Main parts (accessories)**

- PB4200-03172** SHORT RAMP SET PROMEBA
- PB4200-03072** PROMEBA LONG RAMP SET – STD
- PB4200-03310** LONG RAMP SET - 450MM
- PB4200-03271** LONG RAMP SET - 500mm
- PB4200-03290** LONG RAMP SET - 550mm
- PB4205-03090** ULTRA LONG RAMP SET
- PB4100-05150** LONG RAMP SET 900mm

N°	DESCRIPTION OF THE COMPONENTS
1	Stainless steel ramp
2	Roller
3	Ramp inclination regulator

**05.4. Technical specifications / data sheet (accessories)**


Models	PB4200-03172	PB4200-03072	PB4200-03271	PB4200-03271	PB4200-03290	PB4200-03290	PB4100-05150
Length	180mm	350mm	450mm	500mm	550mm	890mm	900mm
Width	581 / 619mm	581 / 619mm	581 / 619mm	581 / 619mm	581 / 619mm	581 / 619mm	581 / 619mm
Roller	NO	Yes	Yes	Yes	Yes	Yes	Yes
Weight	2,86kg	5,35kg	6,9kg	7,4kg	9,4kg	11,5kg	12,5kg

### 06.1. Transport and storage

Before transporting the device, make sure that it is correctly packaged.

#### Transport with a crane (high load)


Thread tapes through the upper platform for its transport (as shown in graphs). Transport the load levelled following all the precepts and regulations for the transport of suspended loads.

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#### Transport with forklift

Place a pallet or flat surface under the rail and transport transversally. Do not load the rail directly on the forklift without a pallet or flat support surface. Do not load the rail longitudinally on the forklift if you do not have a sufficiently long base. Each rail has been thoroughly inspected upon leaving the factory. To ensure that it has not been damaged during transport, please carefully inspect the exterior and interior, and if you find any damage, immediately notify the Promeba or the authorized distributor. For a correct transport operation, the rail box must be levelled.

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To unpack the device, place the box on a flat, stable surface and carefully remove the packing tape. Remove the device from the inside of the case carefully to prevent shocks to the device. Please retain the original packaging for use in case of any other transportation and for storage.

Product damage caused during shipping and handling is not covered by the warranty. Repair or replacement of damaged parts is the responsibility of the customer.

The device should be stored in a cool, dry place, away from direct sunlight. It must not be in contact with substances or chemical agents that may cause damage and reduce its safety performance.

The storage temperature must be between -20°C to 40°C.

The product must not be manipulated in any way (modification, adjustment, addition, substitution). In such cases, all responsibility for any malfunction or injury caused by the product itself will be disclaimed; furthermore, CE certification and product warranty will be considered void.

### 06.2. Preparation

When the product is received:

- The packaging is intact and has no signs of impacts, falls and isn't wet.
- Remove the packaging and expose the material so that all components are visible.
- Check that all the components/pieces on the section 5 are present.

**06 OPERATION**

REV. 2024/12

Before the installation inside the ambulance:


- Ensure that the vehicle in which you want to install the device, is manufactured in compliance with EN 1789.
- Check that the ambulance is capable to store the stretcher support, allowing its movements in all directions in case the device incorporates movements.
- The installation floor is perfectly flat. Non-flat surfaces or parts of the floor that are not perfectly levelled could compromise the correct operation of the product.
- The distance between the ambulance door and the end of the rail is respected. Failure to respect this distance could interfere with the movements of the device with consequences for your safety.

Below are details of each of the anchor points where the rail must be fixed to the floor of the ambulance:


ANCHORING RAIL PB4200-00000, (8 anchor points)

ANCHORING RAIL PB4201-00000, (8 anchor points)

ANCHORING RAIL PB4206-00000, (8 anchor points)


ANCHORING RAIL PB4202-00000, (10 anchor points)


## 06 OPERATION

REV. 2024/12

ANCHORING RAIL PB4203-00000, (10 anchor points)

ANCHORING RAIL PB4204-00000, (10 anchor points)


Before putting the bench into operation, ensure that each anchor point is placed and clasped.

### 06.3. Functioning

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Operate the handle located on the rail right rear side, applying a little effort on it with the foot and accompanying this movement with its slide with the hands.

The total lateral displacement is 300 mm.


### LOADING

- Rail position in active service inside the ambulance (A).
- Open the loading ramp (B).
- Load the stretcher being especially cautious not to hit with the sides of the ramp.
- Fold the loading ramp (C).

### UNLOADING

- Open the loading ramp (B)
- Release the stretcher closure by operating the lever (in case of being installed)
- Pull the stretcher. It will slide outward, hold it with the hands.
- Fold the loading ramp (C).


**06 OPERATION**

REV. 2024/12

**SCOOP STRETCHER FASTENING**


**FIXING**

To place the scoop stretcher, slide it upside down, between the anchoring plate of the stretcher and the tray until it stops against the clamping clips. Once the contact is made, firmly tighten the stretcher until it is anchored.

**EXTRACTING**


To remove the scoop stretcher, pull it firmly until it disengages from the retaining clips. For your comfort, we recommend that you hold the stretcher by the roller to do more force.

**SPINE BOARD FASTENING (MODELS PB4204-00000 Y PB4206-00000)**

**FIXING**

To place the scoop stretcher, slide it upside down, between the anchoring plate of the stretcher and the tray until it stops against the clamping clips. Once contact is made, firmly tighten the stretcher until it is anchored.

To store the spinal board, it must be placed in the lower housing and fixed using the velcro straps for this purpose.

**EXTRACTING**


To remove the scoop stretcher, pull it firmly until it disengages from the retaining clips. For your comfort, we recommend that you hold the stretcher by the roller to do more force.

You can remove the spinal board by undoing the velcro fastening straps.

**TENSORS (PB4200-00000, PB4201-00000 Y PB4206-00000)**

To adjust the tension of the drive cables, the rail has two sets of tensioners for regulation.

Double set of tensioners located on each of the floor fixing plates. Tighten or untighten the cables by screwing or unscrewing the tensors

Manual operation by pedal.

## **07 GENERAL MAINTENANCE**    *REV. 2024/12*

### **07.1.    *Cleaning***

It is essential to keep the equipment clean to ensure proper use and durability of the set. Thorough cleaning should be carried out periodically, especially in areas exposed to dirt that can accumulate dirt such as corners or holes.

Do not use high pressure cleaning systems or bleach to disinfect the product, they can damage it. Instead, use a hydroalcoholic-based disinfectant, wash it with water and let it dry naturally, do not use direct heat sources to dry.

Non-compliance with cleaning operations can lead to the risk of cross infections due to the presence of secretions and/or residues. During all control and disinfection operations, the operator must wear appropriate personal protective equipment, such as gloves, glasses, etc.

When the cleaning procedure includes a disinfectant product, carefully follow the instructions of the manufacturer of the product used regarding the method of application and contact time. Make sure that you have taken all the necessary precautions to ensure that there are no risks of cross-infection or contamination of patients and operators.

Failure to carry out the correct cleaning routine could increase the risk of cross infection, due to presence of body fluids and/or residuals.  
The operator must always wear adequate personal protection such as gloves and mask, etc. during all checking and cleaning procedures.

### **07.2.    *Lubricating***

In general, all moving parts must be greased. Our products leave the factory completely greased and lubricated, however, it is possible that over time and use of the product the elements become degreased, either due to loss of lubrication or due to dirt.

Periodically clean and oil affected areas according to manufacturer's specifications. Check for loose, missing or worn parts.

Periodically inspect all moving parts to ensure that components are tight.

### **07.3.    *Wear areas***

Periodic inspection of the different components of the system for signs of wear is a preventative measure that can reduce breakdowns. Check especially the moving parts such as the ramp and the lateral shift pedal, if present.

In the event that a worn component needs to be changed, you should contact our sales department for more information on ordering spare parts and their installation.

The repair or replacement of product components must necessarily be carried out by qualified personnel, that is, who have received training from the Manufacturer, Promeba S.L. and original spare parts should always be used.

Promeba S.L. declines any responsibility for any damage, direct or indirect, resulting from improper use of spare parts and/or any repair intervention carried out by unauthorized persons.

**07 GENERAL MAINTENANCE**    *REV. 2024/12*

**07.4. Maintenance program**

Although the frequency of the controls is determined by factors such as legal requirements, the type of use, the frequency of use, the environmental conditions during use and storage, it is placed a recommendation on the maintenance that the equipment must have for correct operation throughout its useful life.

All maintenance and revision activities must be recorded and documented with the corresponding technical intervention reports, and the documentation must be kept for at least 10 years from the end of the product's useful life and must be made available to the competent Authorities and/or from the Manufacturer, when requested.

<b>Maintenance review</b>	<b>Every use</b>	<b>When necessary</b>	<b>Every month</b>	<b>Every year</b>
Disinfect	X		X	X
Clean		X	X	X
Inspect	X	X	X	X
Grease the ramp and the lateral displacement				X
Change the scoop stretcher clips				X
Check the tension of the lateral displacement cables		X		

For both the disinfection and cleaning procedures, point 7.1 of this manual must be followed.

During the inspection procedure, which must be done after each use, the following tasks must be carried out:

- Check that all components are present.
- Check the integrity of the device, that is, that there are no breaks, cracks, holes or cuts.
- Check the wear of the product following section 7.3 of this user guide.
- Check that the moving parts slide correctly.
- Check that the accessories provided are present, intact and that they function correctly.

**08 SPARE PARTS**


REV. 2024/12

To request a spare part, indicate:

- Number and description of the set in which the necessary part is found (see table)
- Point out the piece in the drawing, or the area where it is located.

Following these instructions, our sales department will provide you with a detailed drawing of the assembly so that you can specify the reference of the required spare part.

PB4200-00000


**08 SPARE PARTS**

REV. 2024/12

Nº	COMPONENTS DESCRIPTION
1.01	SUPLEMENT SUBASSEMBLY
1.02	PLATFORM SUBASSEMBLY
1.03	RAMP TURNING BUSHING
1.04	STAINLESS STEEL TRAY
1.05	RAMP STOP
1.06	SCOOP STRETCHER CLIP
1.07	LATERAL DISPLACEMENT HANDLE
1.08	TENSOR STOP SUBASSEMBLY
1.09	SCOOP STRETCHER STOP
1.10	SCOOP STRETCHER SUPPORT BUSHING
1.11	SCOOP STRETCHER U STOP
1.12	DAMPER SUPPORT
1.13	PB-417 RAMP SUPPORT
1.14	GAS DAMPER
1.15	HANDLE RUBBER
1.16	BELT FIXATION SUBASSEMBLY
1.17	14x4 HANDRAILSCOOP

Nº	COMPONENTS DESCRIPTION
1.18	SCOOP STRETCHER CLAMPING TAPE
1.19	LARGE ROLL SUBASSEMBLY
1.20	ANCHORAGE SUBASSEMBLY
1.21	TENSOR CABLE SUPPORT
1.22	TENSOR
1.23	CLAMP FIXATION BOLT
1.24	FIXATION BOLT
1.25	D1481 4x18 BOLT
1.26	COMPRESSION SPRING
1.27	BALLS BEARING
1.28	LATERAL DISPLACEMENT GUIDE
1.29	PLASTIC STOPPER 100x60
1.30	PLASTIC STOPPER 80x30
1.31	MUELLA ESPIRAL HILO PIE
1.32	RAMP
1.33	ROLLER BUSHING RAMP
1.34	PLASTIC STOPPER 30x30

PB4201-00000


Nº	COMPONENTS DESCRIPTION
2.01	PLATFORM SUBASSEMBLY
2.02	RAMP
2.03	STAINLESS STEEL TRAY
2.04	DAMPER SUPPORT
2.05	RAMP SUPPORT
2.06	GAS DAMPER
2.07	LARGE ROLL SUBASSEMBLY
2.08	PLASTIC STOPPER 80x30
2.09	ROLLER BUSHING RAMP
2.10	ANCHORAGE SUBASSEMBLY
2.11	TENSOR CABLE SUPPORT
2.12	TENSOR

Nº	COMPONENTS DESCRIPTION
2.13	CLAMP FIXATION BOLT
2.14	FIXATION BOLT
2.15	D1481 4x18 BOLT
2.16	COMPRESSION SPRING
2.17	BALLS BEARING
2.18	LATERAL DISPLACEMENT GUIDE
2.19	RAMP TURNING BUSHING
2.20	PLASTIC STOPPER 30x30
2.21	MUELLE ESPIRAL HILO PIE
2.22	LATERAL DISPLACEMENT HANDLE
2.23	TENSOR STOP SUBASSEMBLY
2.24	HANDLE RUBBER


PB4203-00000

Nº	COMPONENTS DESCRIPTION
4.01	PLATFORM SUBASSEMBLY BÁSICA
4.02	RAMP
4.03	STAINLESS STEEL TRAY
4.04	DAMPER SUPPORT
4.05	RAMP SUPPORT
4.06	GAS DAMPER

Nº	COMPONENTS DESCRIPTION
4.07	LARGE ROLL SUBASSEMBLY
4.08	PLASTIC STOPPER 80x30
4.09	ROLLER BUSHING RAMP
4.10	PLASTIC STOPPER ø25
4.11	RAMP TURNING BUSHING

PB4204-00000

N°	COMPONENTS DESCRIPTION
5.01	SUPLEMENT SUBASSEMBLY
5.02	RAMP TURNING BUSHING
5.03	STAINLESS STEEL TRAY
5.04	RAMP STOP
5.05	SCOOP STRETCHER CLIP
5.06	SCOOP STRETCHER STOP
5.07	SCOOP STRETCHER SUPPORT BUSHING
5.08	SCOOP STRETCHER U STOP
5.09	DAMPER SUPPORT
5.10	PB-417 RAMP SUPPORT

N°	COMPONENTS DESCRIPTION
5.11	GAS DAMPER
5.12	BELT FIXATION SUBASSEMBLY
5.13	14x4 HANDRAIL
5.14	SCOOP STRETCHER CLAMPING TAPE
5.15	LARGE ROLL SUBASSEMBLY
5.16	PLASTIC STOPPER 100x60
5.17	PLASTIC STOPPER 80x30
5.18	RAMP
5.19	ROLLER BUSHING RAMP

## **11 LEGAL NOTICES**

*REV. 2024/12*

- This document may contain technical inaccuracies or typographical mistakes.
- Changes are periodically added to the information herein; these changes will be incorporated in new editions of the publication.
- Promeba, S.L. reserves the right to make any modification or improvement in the products described in this publication if it is appropriate.
- Promeba, S.L. may have patents or patent applications which address themes described in this document. The possession of this document does not entitle no license to these patents. The information contained in this document does not affect or change the specifications or product specifications or warranties of Promeba, S.L.
- No part of this document shall operate as express or implied license or compensation under the intellectual property rights of Promeba, S.L. or .third party.
- All the information contained in this document has been obtained in specific environments and it is presented as an illustration. The results obtained in other operating environments may vary.
- Promeba, S.L. can use or distribute the information that the customer provides in whatever way they see fit, without incurring any obligation with the client.

## 12 PRODUCT WARRANTY REV. 2024/12

Promeba, S.L. guarantees that its products have satisfactorily passed all established quality controls, both functional and material. The duration of the guarantee is 2 years from the date of purchase of the product.

This guarantee will be valid only when it is presented with the original invoice or proof of purchase (indicating the date of purchase, model and the name of the distributor) together with the defective product during the period covered by the guarantee. Promeba, S.L. reserves the right not to offer the free warranty service if the indicated documents are not presented or if the information they contain is incomplete or illegible.

1. This warranty will not apply if the model name or serial number of the product has been altered, erased, disappeared or becomes illegible.

2. This guarantee does not cover transportation costs or risks derived from the transportation of your product to and from Promeba, S.L.

3. This warranty does not cover any of the following:

- a) Periodic maintenance and repair or replacement of parts derived from normal wear and tear.
- b) Expendable material (components that are expected to need periodic replacement during the life of the product, such as belts, straps, etc.).
- c) Damages or defects derived from the improper use, operation or treatment of the product and not due to normal use of the product.
- d) Damages derived from:
  - i. Misuse, including:
    - Treatment that results in damages or physical, superficial or appearance changes of the product.
    - Installation, use or storage of the product in a way that does not

respect the instructions described by Promeba, S.L.

- Maintenance of the product in a way that does not respect the instructions of Promeba, S.L. for proper maintenance.

- Installation or use of the product in a way that does not respect the technical or safety regulations of the country where it is used or installed.

ii. Use of components not provided with the product or incorrect installation of accessory parts not previously tested.

iii. States or defects of the system in which the product is used or incorporated with the exception of other products Promeba, S.L. designed for use with the product.

iv. Use of the product with accessories, peripheral units and other products of a type, condition or standards not established by Promeba, S.L.

v. The manufacturer or distributor will be solely responsible for deciding whether to send the parts for repair, or to replace the product in its entirety. In no case will operators be sent to repair or replacement of the product.

Except in the cases mentioned above, Promeba, S.L. will make no warranties in relation to product, performance, accuracy, reliability, or fitness for logic purpose of the equipment or other type. If this exception is not legal or contemplated by current law, Promeba, S.L. will limit or exclude your warranties only to the extent permitted by applicable law.

The only obligation on the part of Promeba, S.L. in connection with this warranty is to repair or replace parts subject to the terms and conditions of this warranty.

Promeba, S.L. is not responsible for loss or damage to products, this warranty or others, including economic loss or non-assessable

## **12 PRODUCT WARRANTY** *REV. 2024/12*

damage; the price paid for the product; loss of profits, income, information, usufruct or use of the product or associated products or indirect, accidental or critical loss or damage.

This clause refers to whether the loss or damage is due to deterioration or inoperability of the associated product due to defects or unavailability of Promeba, S.L., which has caused downtime, loss of user time or business interruption.

In cases where the law prohibits or limits these liability exclusions, Promeba, S.L. will exclude or limit his liability only to the extent permitted by applicable law.

For example, there are countries that prohibit the exclusion or limitation of damages caused by negligence, reckless negligence, willful misconduct, fraud and similar actions. The responsibility of Promeba, S.L. In this guarantee will not exceed, in any case, the price paid for the product, but if the current law allows only limitations of greater responsibilities, these will apply.

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*All rights reserved. Variations can be done whitout notice.  
Promeba, S.L. is to be considered not responsible for damages caused by the lack or the wrongness of the information here mentioned.*

SCAN QR CODE

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Ctra C-16 Km 59.5 · 08650 Sallent (Barcelona) · SPAIN

T. 93 837 12 00

[promeba@promeba.com](mailto:promeba@promeba.com) · [www.promeba.com](http://www.promeba.com)