MANUAL DE INSTRUCCIONES / USER GUIDE
SILLA PS-251 + PA-260 / CHAIR PS-251 + PA-260

DT03_02_05 REV. 2026/07

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01	MODELS.....	56
01.1	Products.....	56
01.2	Accessories.....	56
02	INTRODUCTION.....	57
02.1	Using this guide.....	57
02.2	Legend of symbols.....	57
02.3	Service request.....	58
02.4	Life span.....	58
02.5	Demolition.....	59
02.6	Reference standards.....	59
02.7	Labelling.....	60
03	INTENDED USE.....	61
03.1	Intended use.....	61
03.2	Contraindications and side effects.....	61
03.3	Physical requirements of the operators.....	61
04	WARNINGS.....	62
04.1	General warnings.....	62
04.2	Specific warnings.....	64
04.3	Residual risk.....	65
05	DESCRIPTION OF THE PRODUCT.....	66
05.1	Main components.....	66
05.2	Technical data sheet.....	69
05.3	Main components (accessories)	72
05.4	Technical data sheet (accessories)	75
06	OPERATION.....	77
06.1	Transport and storage.....	77
06.2	Preparation.....	77
06.3	Functioning.....	81
06.4	Functioning of the track.....	86
07	GENERAL MAINTENANCE.....	94
07.1	Cleaning.....	94
07.2	Lubrication.....	94
07.3	Wear areas	94
07.4	Maintenance program.....	95
08	SPARE PARTS.....	96
09	TRAINING REGISTER.....	101
10	MAINTENANCE REGISTER.....	102
11	LEGAL NOTICES.....	103
12	PRODUCT WARRANTY.....	104

01	MODÈLES.....	56
01.1	Produits.....	56
01.2	Accessoires.....	56
02	INTRODUCTION.....	57
02.1	Utilisation du manuel.....	57
02.2	Légende des symboles.....	57
02.3	Service client.....	58
02.4	Durée de vie.....	58
02.5	Élimination.....	59
02.6	Réglementation applicable.....	59
02.7	Étiquetage.....	60
03	INTENDED USE.....	61
03.1	Intended use.....	61
03.2	Contraindications and side effects.....	61
03.3	Physical requirements of the operators.....	61
04	WARNINGS.....	62
04.1	General warnings.....	62
04.2	Specific warnings.....	64
04.3	Residual risk.....	65
05	DESCRIPTION OF THE PRODUCT.....	66
05.1	Main components.....	66
05.2	Technical data sheet.....	69
05.3	Main components (accessories).....	72
05.4	Technical data sheet (accessories).....	75
06	OPERATION.....	77
06.1	Transport and storage.....	77
06.2	Preparation.....	77
06.3	Functioning.....	81
06.4	Functioning of the track.....	86
07	GENERAL MAINTENANCE.....	94
07.1	Cleaning.....	94
07.2	Lubrication.....	94
07.3	Wear areas.....	94
07.4	Maintenance program.....	95
08	SPARE PARTS.....	96
09	TRAINING REGISTER.....	101
10	MAINTENANCE REGISTER.....	102
11	LEGAL NOTICES.....	103
12	PRODUCT WARRANTY.....	104

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REV. 2026/07

01.1. Products

REFERENCE	UDI CODE	TRADE NAME
PS2510-00000	8437001519125	EVACUATION CHAIR PS-251
PA2600-00000	8437001519163	ELECTRIC SYSTEM TRACK PA-260
PA2650-00000	8437001519170	MANUAL TRACK PA-265
PS2510- PA2600	8437021314670	ELECTRIC CHAIR PS-251 + PA260

01.2. Accessories

REFERENCE	UDI CODE	TRADE NAME
PA2810-00000	8437001519224	LOCK SET PA-281
PA2710-00000	8437001519231	LOCK SET PA-271
PA5770-00000	8437001519200	LOCK SET PA-577

02 INTRODUCTION

REV. 2026/07

02.1. Using this guide

The manual provides instructions for use and maintenance of the product, as well as technical aspects, operation, spare parts and safety.

It is recommended before the operation of the product to read carefully this manual in order to avoid damages caused by a misuse, regardless of the level of experience acquired in the past with similar devices.

Do not lose this document. It should be accessible to any doubt that could appear by medical personnel. In case the device has some instructions for use of a different product than the one received, it is necessary to contact the manufacturer immediately before using the product.

Remember that a good use and maintenance are necessary for the proper operation of the product.

Each product includes an identifying adhesive label with the serial number and model. Keep these numbers so that you can refer them to the dealer if necessary.

02.2. Legend of symbols

SYMBOL	EXPLANATION / DESCRIPTION
	MANUFACTURER symbol. This symbol is accompanied by the name and address of the manufacturer, adjacent to the symbol (PRODUCTOS METÁLICOS DEL BAGÉS S.L., Ctra. C-16 Km 59.5, 08650 Sallent (Barcelona)).
	Indicates the manufacturer's reference number to identify the medical device. PROMEBA, S.L. uses this symbol to set each internal reference for each configuration and business variant
	Indicates the manufacturer's serial number to identify a specific medical device.
	Indicates the manufacturing date. The symbol must be accompanied by a manufacturing date (yyyy-mm), adjacent to the symbol.
	It is placed to inform that the product is a "Medical Device"
	CE symbol without the intervention of a Notified Organism, as a Medical Device classified as Class I according to the EU Regulation 2017/745 on Medical Devices.
	Symbol for the Unique Device Identifier.
	Symbol "See instructions for use or operating instructions"
	Symbol "caution". This symbol is placed to warn of the need for the user to refer to important precautionary information in the operating instructions, such as warnings and cautions not otherwise found on the label.
	Waste symbol for electrical and electronic equipment "WEEE" (Waste Electrical and Electronic Equipment). Recycle: Electronic equipment. DO NOT THROW GARBAGE
	"Dangerous voltage" symbol. To warn of electricity

02 INTRODUCTION

REV. 2026/07


Symbol "Caution". For a general warning.


Warning, crushing of hands

"Double insulation" symbol. The protection of appliances marked with this symbol is ensured by double insulation and does not require a safety electrical connection to earth (ground).

Cot duty cycle: 16.7% (less than 60 seconds on, more than 300 seconds off)

Indicates the temperature limits. The upper and lower temperature limits should be indicated adjacent to the horizontal lines.

02.3. Service request

For information of the correct interpretation of the instruction manual, the use, maintenance, installation and restoration of the product, please contact Promebea customer service: T. 93 837 12 00, email: promebea@promebea.com or write to PROMEBEBA, S.L. - Ctra C-16 Km 59.5 · 08650 Sallent (Barcelona) · SPAIN.

Whenever you want to communicate any matter related to a product, you will be asked for the product reference (**REF**) and the serial number (**SN**). It is strongly recommended to keep the label containing this information and found on the original packaging of the product.

02.4. Life span

If used as described in the following instructions, this device has a useful life of 10 years from the date of purchase.

This useful life can be extended with annual reviews carried out by the manufacturer, which uses specialized and authorized internal and external technicians.

In case these annual checks are not carried out, the device must be disposed of according to the information in paragraph 02.5 and the manufacturer must be notified.

Only the manufacturer or an authorized center can extend the life of the device, if it meets the safety and legal requirements.

Promebea, S.L. will not accept any responsibility for malfunction or damage caused by the use of devices that have not been checked by the manufacturer or authorized center, or that have exceeded the maximum allowed useful life.

02 INTRODUCTION

REV. 2026/07

02.5. Demolition

When the devices are no longer suitable for use, if they have not been contaminated by any particular agent, they can be disposed of as electric waste, otherwise, follow the current demolition regulations.

WARNINGS FOR THE CORRECT DISPOSAL OF THE PRODUCT IN ACCORDANCE WITH EC DIRECTIVE 2012/19/EU WEEE:

At the end of his life, the product must not be disposed as household waste. Can be taken to special recycling centers provided by local government, or return it to the dealer on purchase of a new device of the same type and used for the same functions. Dispose of the product separately avoids possible negative consequences for the environment and human health resulting from inappropriate disposal and allows to recover the materials in order to obtain significant savings in energy and resources. The symbol on the label indicates separate collection of electrical and electronic equipment.

02.6. Reference standards

REFERENCE	TITLE OF DOCUMENT
EU Regulation 2017/745	EU Regulation on Medical Devices
EN 1789:2007+A2:2014	Medical vehicles and their equipment - Road ambulances
UNE-EN 1865-4:2012	Patient handling equipment used in road ambulances - Part 4: Foldable patient transfer chair
EN 60601-1:2006/A1:2013 IEC 60601-1:2005/A1:2012	Medical electrical equipment - Part 1: General requirements for basic safety and essential performance
EN 61000-6-1:2007 IEC 61000-6-1:2005	Electromagnetic compatibility (EMC) - Part 6-1: Generic standards - Immunity for residential, commercial and light-industrial environments
EN 61000-6-3:2007+A1:2012 IEC 61000-6-3:2006+A1:2010	Electromagnetic compatibility (EMC) -- Part 6-3: Generic standards - Emission standard for residential, commercial and light-industrial environments


As a distributor or end user of the products manufactured and/or distributed by Promeba, S.L., you are strictly obliged to know the legal provisions in force in the country of destination of the goods, applicable to the devices to be supplied (including the regulations relating to technical specifications and/or safety requirements) and, therefore, understand the necessary requirements to ensure compliance of the products themselves with all the legal requirements of the territory.

02 INTRODUCTION

REV. 2026/07

02.7. Labelling

Each product incorporates an identification label, placed on the device itself and/or on the original packaging of the product. It must never be removed or covered. This label includes the data of the company which produced the device () the product code / reference () of the product, the manufacturing date () the serial number () the UDI-PI code () and the CE distinctive (). Please keep these numbers so you can inform the dealer if necessary.

The UDI-PI code is made up of the following fields:

- (01)08437001519163: UDI code of the product.
- (11)231127: it gives information about the batch release date of the product following the format *aa**mm**dd*, where *aa* is the batch release year, *mm* is the month of the batch release and *dd* is the day of the batch release.
- (21)01554823001: identifies the serial number of the product giving information about the manufacturing order (0155) including the week when the product was produced (4823) and a correlative number depending on the units that has been produced in that production order (001).

03 INTENDED USE

REV. 2026/07

03.1. *Intended use*

The product ELECTRIC CHAIR is a device indicated to effortlessly transport patients (over 40kg) in hospital and pre-hospital environments thanks to its electric transport function when going up and down stairs.

03.2. *Contraindications and side effects*

Any contraindication or side effect have been resulted from the clinical analysis of the ELECTRIC CHAIR family of products.

03.3. *Physical requirements of the operators*

Promeba electric chairs are intended exclusively for professional use. Each operator must be trained to transport patients safely and efficiently. Do not allow untrained persons to assist during use of the product, as they may cause injury to others or themselves. The device is intended to be used by medical transport technical personnel or ambulance technicians who have knowledge and experience with patient transport.

Operators must be able to provide the necessary assistance to the patient.

Users must be able to safely lift and manipulate the weight of the assembly consisting of the chair, the patient, and any other equipment used with the device.

For patient loading techniques, for particularly heavy patients, for interventions on steep terrain or in particular and unusual circumstances, the presence of several operators is recommended (not just 1 as expected in standard conditions).

Operators must be able to evaluate any anomaly of the system/assembly, communicating the problem to those responsible and interrupting the use of the devices if necessary to ensure patient safety.

The capacities of the various operators must be considered before determining his role in the employment of the board.

04 WARNINGS

REV. 2026/07

04.1. General warnings

1. The product must be used by trained personnel only, having attended specific training for this device and not for similar products.
2. Establish a maintenance and periodic control program and identify a reference manager. The person in charge of the ordinary maintenance of the device must guarantee the basic requirements provided by the manufacturer in these instructions for use. It is recommended that the subject check at least every 6 months if there are updated user manuals available or if there are any changes related to his / her product. This information is freely available on the website <http://promeba.com/>.
3. Training routines must be registered on a special register in which the names of those trained, of the trainers, date and place are indicated. This register, which will certify the eligibility of the operators to use the Promeba, S.L. device, has to be kept for a period of 10 years after the disposal of the device itself. This register will be made available to the competent authorities and/or manufacturer if requested.
4. Promeba, S.L. is always at your disposal to plan trainings on products.
5. Before carrying out any kind of operation on the appliance (training, installation, use), the operator must carefully read the enclosed instructions, paying particular attention to the correct safety precautions and to the procedures to be followed for installation and for correct use.
6. If the instructions belong to another device and not the device received, inform the manufacturer immediately and avoid use of the device.
7. In the case of any doubts as to the correct interpretation of the instructions, please contact Promeba, S.L. for any necessary clarifications.
8. Do not allow untrained persons to help during the use of the device, because they could cause damage to the patient or to themselves.
9. Before each use of device the perfect operating state of the device must be checked as specified in the instruction manual. If any damage or abnormalities which could in any way influence the correct functioning and the safety of the device, of the patient and/or of the user are detected, the device must be immediately removed from service and the manufacturer must be contacted.
10. If any failure or incorrect functioning of the device is detected, it must be immediately substituted with a similar item so that the rescue procedures are guaranteed without any interruption.
11. Use of the device in anyway other than described in this manual is forbidden.
12. Do not alter or modify in any way the device; any such interference could cause malfunctions and injury to the patient and/or rescuer.
13. The device must not in any way be tampered with (modification, adjustment, addition, replacement). In such cases all responsibility will be denied for any malfunctions or injuries caused by the device itself; moreover CE certification and product warranty will be considered void.

04 WARNINGS

REV. 2026/07

14. Use only original or PROMEBA approved components/spare parts and/or accessories to carry out all operations without causing alterations or modifications to the device; Otherwise, any responsibility for the operation is declined.
15. Ensure that all the necessary precautions are taken in order to avoid the hazards that can arise as the result of contact with blood or body fluids.
16. Register and store with these instructions: lot number, place and date of purchase, first date of use, date of checks, name of users and any other comments.
17. When the device is being used, the assistance of qualified staff must be guaranteed.
18. The product must not be exposed or come into contact with thermal sources of combustion and flammable agents, but must be stored in a dry, cool place, protected from light and sunlight.
19. Store and transport device in its original packaging.
20. Do not store the product under other more or less heavy materials, which can damage the structure of the product.
21. Attention: laboratory testing, post production tests and instruction manuals cannot always consider every possible scenario for use. This means that in some cases the performance of the product could be notable different from results to date obtained. Instructions are continually being updated and are under tight surveillance of fully qualified staffs with adequate technical formation.
22. As a distributor or end users of products manufactured and/or marketed by Promeba, S.L., you are strictly required to have a basic knowledge of any legal requirements existing in the goods final destination country applying to the devices contained in this supply (including laws and norms regarding technical specifications and/or safety requirements) and therefore you are also strictly required to have the necessary knowledge to guarantee all aspects regarding the total conformity of the products to the regulations in the relevant territory.
23. Promptly notify PROMEBA regarding any revisions to be made by manufacturer in order to guarantee the conformity of the product to the territory's legal specifications (including those resulting from rules and/or norms of other nature).
24. Act with all due care and diligence, and contribute to ensure conformity to general safety requirements of all devices marketed in the territory, by providing final users with all necessary information for carrying out periodical checks on their devices, as specified in the relevant user manual.
25. Actively contribute to product safety checks on products sold, by communicating any relevant risk analysis information both to the manufacturer and to any competent authorities so that the necessary actions can be promptly taken.
26. Be aware that in the event of any failure to conform to the above mentioned requirements you will be deemed fully responsible for all damages that might occur. Therefore we expressly disclaim any responsibility and/or liability for your non-compliance with the present "Regulatory provisions".
27. Follow internal procedures and protocols approved by your organization.

04 WARNINGS

REV. 2026/07

28. Follow approved Emergency Medical Service procedures for patient immobilization and transportation.

04.2. Specific warnings

1. Use only accessories/spare parts that are original or approved by Promeba, S.L. when performing any operation so as not to cause any alteration or modification to the device, otherwise we do not assume any responsibility for the correct operation or resulting damage of the device to the patient or the operator and the guarantee will be considered void, as well as compliance with the Regulation 2017/745.
2. Always respect the maximum capacity of the device, as indicated in this user's manual. Maximum load capacity means the total weight distributed according to the human anatomy. In determining the load of the total weight on the product, the operator must consider the weight of the patient, the equipment and the accessories. Moreover, the operator must consider that the overall dimensions of the patient do not reduce the functionality of the device.
3. Inside the ambulance, place and adjust the device taking care not to obstruct other parts or other equipment. It is recommended not to open the equipment inside the ambulance, but rather to do it outside the vehicle.
4. Never leave the patient unassisted on the device, because he may be injured.
5. Do not use bleach to disinfect the device. Use a hydroalcoholic based disinfectant and wash with water.
6. The device and all its components, after washing, should be allowed to dry completely before storing.
7. Avoid contact with sharp objects.
8. Do not use the device if it is pierced, deformed or shows signs of loss of integrity.
9. Do not operate in case the weight has not been distributed correctly.
10. Condensation, water, ice and dust accumulations can affect the correct operation of the device, making it unpredictable.
11. Always immobilize the patient at least with the belts supplied by the manufacturer during any movement of the stretcher, including loading and unloading to the ambulance.
12. Do not use the device as a support for climbing or leaning on and reaching high points in the ambulance. The product may be damaged and the technician or patient may be harmed.
13. Disinfection activities must be carried out in accordance with the validated cycle parameters, indicated in the specific technical standards.
14. Promeba S.L. shall not be liable if unauthorized / incompatible stretcher fastener is used.
15. Do not leave the batteries uncharged for more than 2 MONTHS. They could break down.

04 WARNINGS

REV. 2026/07

16. Incorrect loading of the chair can cause damage and/or injuries. Use only batteries approved and authorized by Promeba S.L. and recharge only with chargers approved and authorized by Promeba S.L.
17. Do not use the chair with the track carriage without being properly locked as this may cause damage and/or injury. Make sure you have properly locked the track carriage before moving the patient.
18. Do not exceed the maximum capacity allowed by the chair. Exceeding it may cause deformations to the product or harm to the patient.
19. Do not lubricate the chair when it is connected to the power.
20. In repeated up and down maneuvers of up to 10 uninterrupted cycles, a mandatory stop must be made for at least 15 seconds so that the electric chair's motor-reduction system recovers.
21. Always use both hands to transport the chair to avoid the risk of tipping over.
22. Pay full attention to the team and the patient during the performance of the service.

04.3. *Residual risk*

The residual risks listed below have been identified only with reference to the intended use of the device:

1. A patient without assistance is at risk of injury. Never leave the patient alone and always help them while using the chair.
2. The operation carried out by unqualified personnel could cause harm to both the patient and the personnel performing the operation.

05.1. Main components

PS2510-00000 EVACUATION CHAIR PS-251

N°	DESCRIPTION OF THE COMPONENTS	N°	DESCRIPTION OF THE COMPONENTS
1	Foldable Armrests	5	Front swivel wheels with brake
2	Leg and chest belts with adjustable safety lock	6	Rear wheels with brake
3	Telescopic handles with the possibility of total extraction	7	Foldable transport handles
4	Integrated folding footrest	8	Track anchor system

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PA2650-00000 MANUAL TRACK PA-265

N°	DESCRIPTION OF THE COMPONENTS	N°	DESCRIPTION OF THE COMPONENTS
1	Handle coated with non-slip cover	5	Inclination regulation system
2	Telescopic handle adjustable system	6	Tensioners
3	Locking system between chair and track	7	Anti friction tracks
4	Tubular structure in painted aluminium	8	Anti roll-over wheels

05.2. Technical specifications / data sheet

Models	PS-251	PA-260	PA-265
Length	700mm	170mm	685/1180mm
Width	520mm	350mm	550mm
Height	1025mm	1090-1520mm	1090mm
Height folded	225mm		
Weight	13kg*	19kg**	9kg
Max load	250kg	250kg	250kg
Color	White	Black	Black
Øfront wheel	120mm	2 x Ø50mm	2 x Ø50mm
Ørear wheel	200mm		
Handles	4 (2del. + 2 tras.)		
Max. Autonomy		140min /130kg	
Voltage		12V	
Consumed potence		30 A	
Loading time		60 min	


*Without belts, armrests, footrests and handles

**without battery

05 PRODUCT DESCRIPTION


REV. 2026/07

PS2510-00000

PA2600-00000


PA2650-00000

05.3. Main parts (accessories)


PA2810-00000 LOCK SET PA-281

N°	DESCRIPTION OF THE COMPONENTS
1	Quick and easy locking system
2	Centering platform
3	Fixing suitable for any flat surface
4	Belt fastening system



PA2710-00000 LOCK SET PA-271

N°	DESCRIPTION OF THE COMPONENTS
1	Quick and easy locking system
2	Clamp for quick anchorage
3	Integrated charging system
4	Support for loading the track
5	Tubular aluminium structure
6	Centering platform

PA5770-00000 LOCK SET PA-577


N°	DESCRIPTION OF THE COMPONENTS
1	Belt fastening system
2	Centering platform
3	Wheels slots

05.4. Technical specifications / data sheet (accessories)

Models	PA-281	PA-271	PA-577
Length	300mm	535mm	370mm
Width	395mm	310mm	185mm
Height	468mm	1157mm	766mm
Weight	11kg	19,8kg	6kg
Suitable for	PS-251 PS251+PA-260	PS-251 PA-260 PS251+PA-260	PS-251 PS251+PA-260
10G	YES	YES	NO


PA2810-00000

05 PRODUCT DESCRIPTION

REV. 2026/07

PA2710-00000

06 OPERATION

REV. 2026/07

06.1. Transport and storage

Before transporting the appliance, make sure that it is correctly packaged. During transport always fix the load. If piling up is necessary, always follow the scheme shown on figure 1. Transport the load levelled and following all precepts and rules for the load transportation, ensuring also that there are no risks of shocks, bumps or falls during the transport itself.

To unpack the device place the box on a flat, stable surface and carefully open the seal. Remove the device from the inside of the box following the scheme shown on figure 2.

Keep the original packaging for use in case of any other transport and for storage.

Damage to the appliance caused during transport and handling is not covered by the guarantee. Repairs or replacement of the damaged parts are the responsibility of the client.

The device must be stored in a dry, cool area away from direct sunlight. It must not be placed in contact with any substances or chemical agent which could cause damage and reduce safety characteristics.



The storage temperature must be between -20°C and 60°C.

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SIZE: 1160 x 550 x 305 mm

WEIGHT: 35 Kg.

The product must not be manipulated in any way (modification, adjustment, addition, substitution). In such cases, all responsibility for any malfunction or injury caused by the product itself will be disclaimed; furthermore, CE certification and product warranty will be considered void.

06.2. Preparation

When the product is received:

- Verify that the packaging is intact and has no signs of impacts, falls and isn't wet.
- Remove the packaging and expose the material so that all components are visible.
- Check that all the components/pieces on the section 5 are present.

The appliance must be checked before every use so as to reveal any working abnormalities and/or damage caused by transport and/or storage. In particular, check:

06 OPERATION

REV. 2026/07

- General functionality of the device.
- Cleaning the device (remember that the lack of cleaning can cause the spread of infections).
- Absence of cuts, holes, breaks in the structure, including the belts.
- Correct fixing of all nuts, bolts and screws.
- Correct attachment of the belts to the chair.
- Correct operation of the belt closure.
- Condition of moving parts, wheels and belts.
- Integrity of the seams.
- There are no tubes or metal sheets that show bends or cracks.
- The backrest and the seat do not show damage or structural cracks.
- The welds are intact, without cracks or breaks.
- Wheels are securely fixed, stable and work properly.
- Wheels are free of dirt or debris.
- Brakes work properly.
- Track belts run and have the correct tension for use.
- Correct operation of the springs.
- Armrests raise and lower correctly.
- The transport handles open, close and lock correctly.
- Telescopic handles open, close and lock properly.
- Presence of all labeling.
- Lubrication of moving parts.
- The emergency vehicle is equipped with the locking system for the PS-190 chair.

If the above conditions are met, the device may be considered ready for use; otherwise it must immediately be removed from service and contact the manufacturer.

Before the installation inside the ambulance:


- Ensure that the vehicle in which you want to install the device, is manufactured in compliance with EN 1789.
- Check that the ambulance is capable to store the chair support.
- The installation floor is perfectly flat. Non-flat surfaces or parts of the floor that are not perfectly levelled could compromise the correct operation of the product.
- The distance between the ambulance door and the end of the support is respected. Failure to respect this distance could interfere with the movements of the device with consequences for your safety.

Below are details of each of the 10 anchor points where the rail must be fixed to the floor of the ambulance:

06 OPERATION

REV. 2026/07


PA2810-00000

06 OPERATION

REV. 2026/07

PA2710-00000

06 OPERATION

REV. 2026/07

06.3. Functioning

ARMRESTS

The armrest has two positions:

1. Lifting for storage (See in Figure 2).
 2. Lowered for use (See in Figure 1).
- To lower the armrest, push it down until it stops.
 - To lift the armrest, lift it until it stops.

Instructions for using the armrest:

- Keep the armrest raised when it is not in use.
- Keep the armrest clean, avoiding patient contact with external fluids.

FIGURE 1 - LOWERED ARMREST

FIGURE 2 - LIFTED ARMREST

BELTS

Use patient restraint systems to help the patient in the chair.

Follow these instructions to properly anchor the belts:

1. Join the two buckles indicated in Figure 3, following the direction of the arrows until the system will be closed and then be sure that it does not open.

Instructions for use of the belts:

- Make sure that they are well secured in their respective anchor points, which can be observed in Figure 4.
- Use a minimum of two belts.
- A safety belt should be secured on the seat of the chair to keep the patient on the chair.

FIGURE 3 - CLOSURE BELTS

FIGURE 4 - BELT POSITION

TELESCOPIC HANDLES

The telescopic handles adjust to two locked positions. Adjust the handles according to the preferences of the operator and the task.

EXTEND THE HANDLES:

1. Press the release button (Figure 5) and pull the handle until you find the locking position.
2. Release the button and slide the handle forward until it clicks into place.
3. Lock both handles in the same position. Verify that both handles are locked when you try to push or remove them without

FIGURE 5 - NON-EXTENDING HANDLE

FIGURE 6 - EXTENDED HANDLE

06 OPERATION

REV. 2026/07

pressing the release buttons. When the locks are engaged, the handles will not move.

RETRACT THE HANDLES:

1. Press the release button (Figure 6) and push the handle until you find the blocking position.
2. Release the button and slide the handle backwards until it clicks into place.
3. Lock both handles in the same position. Verify that both handles are locked when you try to push or remove them without pressing the release buttons. When the locks are engaged, the handles will not move.

FOOTREST

The footrest has two positions:

1. Lifting for storage (See in Figure 7).
 2. Lowered for use (See in Figure 8).
- To lower the footrest, lower it down until it stops.
 - To lift the footrest, lift it until it stops.

Instructions for footrest use:

- Use the footrest correctly to ensure that it does not interfere with the patient's operator feet.
- Before transferring the patient to the chair, make sure that the footrest is raised in the storage position. Sit down the patient before lowering the footrest.
- At the time to move the patient from the chair, lift the footrest to the storage position before loosening the belts. This will prevent the patient from trying to stop before the footrest has been stowed out of the way of his feet.
- Keep the repository elevated when it is not in use.

FIGURE 7 - LOWER FOOTREST

FIGURE 8 - LIFTED FOOTREST

FRONT SWIVEL WHEELS BRAKE

The front wheels of the chair are equipped with brakes to prevent the chair from moving during the transfer of the patient or during a stop.

- Be with the chair and keep control of it at all times. Do not use wheel locks as a substitute for operator control.
- To activate a wheel lock, press down the rear end of the lock pedal. (See in Figure 9).
- To deactivate a wheel lock, press down the front end of the lock pedal. (See in Figure 9)

FIGURE 9 - ROTATING WHEEL BRAKES

06 OPERATION

REV. 2026/07

REAR WHEELS WITH BRAKE

The rear wheels of the chair are equipped with brakes to prevent the chair from moving during the transfer of the patient or during a stop.

- Be with the chair and keep control of it at all times. Do not use wheel locks as a substitute for operator control.
- To activate a wheel lock, press down the rear end of the lock pedal. (See in Figure 10).
- To deactivate a wheel lock, press down the front end of the lock pedal. (See in Figure 10)


FIGURE 10 - REAR WHEELS WITH BRAKES

TRANSPORT HANDLES

The transport handles have two positions:

1. Lowered for storage (See in Figure 12)
 2. Lifted for use (See in Figure 11)
- To lower the carrying handles, press the release button (Figure 11) and push them down until they stop.
 - To lift the carrying handles, lift them until they stop.

Instructions for use of transport handles:

- Keep the carrying handles down when they are not in use.
- Keep the carrying handles clean, avoiding operator contact with external fluids.

FIGURE 11 - LIFTED TRANSPORTATION HANDLES

FIGURE 12 - TRANSPORT HANDLES LOWERED

06 OPERATION

REV. 2026/07

TRACK ANCHOR SYSTEM

Use the track anchor system when you want to attach a track to the chair to facilitate patient transport by stairs.

INSTRUCTIONS FOR ANCHORING THE TRACK:


First, place the track in an extended position, reclining on the tread bands as shown in figure 1.

With the chair unfolded, bring the lower hooks of the chair closer to the track, so that they are aligned with the track anchors, slightly lifting the chair attach the lower anchors as shown in figure 2.

Then, turn the chair on the power track, with a dry blow, you should hear a CLACK sound and the coupling will have been produced correctly, see image 3.

To end, lift the assembly and close the track so that the assembly is according to the image 4.

For the decoupling process, pull the handles located at the top of the track and it will be released.

During the operation of the coupling and decoupling of the track, the chair shall remain at all times without any patient and on a surface as flat possible. If not, both patient and personnel may result in wounds and may result in damaged important product elements.

06 OPERATION

REV. 2026/07

FOLDING AND UNFOLDING THE CHAIR

To fold the chair locate the unlocking cable on the back of the backrest, with one hand holding the chair and one foot in the lower tube that joins the wheels, pull the cord toward you and bring the chair closer to the floor so that the seat approaches the backrest. When the chair is completely folded, the locking shelves must be in the locked position.

To unfold the chair follow the process in reverse and make sure that the tubes return to their natural locking position.

The front wheels of the chair and all handles should remain folded and locked when the folding / unfolding operation is done. Otherwise, you could damage important elements of the chair.

06 OPERATION

REV. 2026/07


06.4. Functioning (track)

START UP OF THE TRACK

To start using the power track, it is necessary to charge the battery. To do this, locate the external socket with female connector 2.1, located on the back of the track as shown in the picture (A).

If you have the PA-271 support and its installation and you have requested the PA-260 track with integrated load, you can start loading the track as shown in the image (B), place the track on the stand and start charging automatically.

Always charge the battery up with the charger supplied by PROMEBA to its maximum charge level (approx. 60min. are enough).

Only charge the batteries of the track with the charger supplied by PROMEBA, otherwise it could damage the batteries and you will lose their warranty.

External charger for PROMEBA battery.

High performance battery for PROMEBA power track.

06 OPERATION

REV. 2026/07

GENERAL WARNINGS FOR THE TRACK

- Unless 1 operator is enough to manage the track, two trained operators and a third person as an “observer” are recommended.
- Verify that the tracking system is completely open and locked before using it.
- Never lubricate the belts. Lubrication can cause the belts to work unpredictably, which can cause injury to the patient and/or operators.
- Moisture, water, snow, ice or debris on or between track and belts can cause uneven belt operation that causes sudden changes in the weight that operators must bear. Make sure the tracks and belts are clean and dry before using the chair on the stairs.
- Humidity, water, snow, ice or debris on ladders can cause poor balance for operators. To avoid possible injuries, clear the stairs or select an alternative route.

FIGURE 19 - UNLOCKING KNOB FOR DEPLOYMENT OF THE TRACK

OPENING THE TRACK SYSTEM

1. Hold the unlocking knob to unfold the track shown in Figure 19 and firmly pull it in the indicated direction until you see that the deployment adjustment system moves, and push the track in the direction shown in Figure 20 , placing it in the desired position (Position 1, Position 2, Position 3) as it suits you, depending on the conditions you are in when moving the chair up the stairs.
2. Verify that the system has been locked by pushing it as if you wanted to close it. If the blockage is completely engaged, the track system will not close.

FIGURE 20 - OPENING OF THE TRACK SYSTEM

CLOSING THE TRACK SYSTEM

1. Hold the release knob to unfold the track shown in Figure 19 and pull firmly in the indicated direction until you see the deployment regulation system move, close the track in the direction shown in Figure 21 to that it is completely closed.
2. Verify that the system has been locked by pushing it as if you wanted to open it. If the blockage is completely engaged, the track system will not open.

FIGURE 21 - CLOSURE OF THE TRACK SYSTEM

06 OPERATION

REV. 2026/07

OPERATION UP AND DOWN

Starting from the transport position and always with a trained sanitary technician:

UP


- Attach the chair to the track as explained in page 33.
- Extend the telescopic handles of the chair as explained in page 82.
- Fold the backrest handles as explained in page 84.
- Extend the footrest once the patient is sitting as explained in page 83.
- Fasten the patient with the two belts supplied as explained in page 82.
- Fix the patient's head using the immobilizer located on the track (A).
- Set the track to "stand by" by pressing the ON / OFF button on the control panel (B).
- Approach the chair to the base of the ladder and gently recline the track drive above the first steps of the ladder.
- If there is a second operator, he/she must extend the telescopic handles before starting the ride and accompany at all times during the course (2).
- Press the "UP" button to go up the stairs (3).
- Once the section of the stairs is finished, it is very important to fold all the necessary elements during the plane transfer (4).

DOWN

- Approach the chair to the beginning of the stairs.
- In case there is a second operator, he/she should help placing the chair at the beginning of the ladder and should assist the inclination of the track.
- Press the "DOWN" button to go down the stairs (3).

The patient should stay fixed at all times when the lift up / down is operated.
It is not permitted in any case to mobilize patients without complying with all the safety measures. Established.


06 OPERATION

REV. 2026/07

CONTROL PANEL FUNCTIONS

All of the track's electrical functions are controlled from the touch panel located at the top of the track.

See its functions in the following diagram:


CAPTION OF ERRORS IN THE CONTROL PANEL

The LEDs of battery status indicate possible errors of the system. In the following table we describe the error according to its combination:

1. CONNECTION ERROR BETWEEN THE TRACK AND THE SUPPORT
2. ERROR OF SUPPLY ELECTRIC IN THE SUPPORT
3. OVERVOLTAGE ERROR ON BATTERIES
4. LOAD TIME EXCESS ERROR
5. ERROR IN EXCESSIVE TIME RISE MOVEMENT
6. ERROR IN LOWER MOVEMENT FOR EXCESSIVE TIME
7. ERROR IN THE TRACK CONTROL DRIVER
8. ENGINE ERROR WHEN EXCEEDING MAXIMUM INTENSITY
9. LOW VOLTAGE MOTOR ERROR
10. OVERLOAD ERROR


06 OPERATION

REV. 2026/07

ADJUST THE TRACK TO THE STAIRS

The PA-260 track incorporates a manual regulator that allows adjusting the slope of the track according to the slope of the ladder. In this way an ergonomics of use is allowed during the maneuver of up and down, being the position of the patient always of about 15° respect to the floor.

There are 3 settings available:

Pull the latch as shown in the chart on the label and adjust the incline of the track according to the slope of the ladder.

Always make sure that the regulator is securely locked before starting up or down the ladder.

BATTERY

The battery charging time is 5-6 hours. When the battery is charging, the charger light is red. When fully charged, the light turns green. It is recommended to continue charging for an additional hour after the green light has come on.

The battery can last 1.5 hours once fully charged, which is equivalent to going up and down a 5-story building 30 times. When the chair is not in use, close the battery switch.

If the chair is not used for a long time, the battery should be charged every 3 months. The battery can be charged 500 times. The service life of the battery is 3-5 years.


06 OPERATION

REV. 2026/07

PLATE CIRCUIT BOARD

The below diagram shows the main PCB of the track PA-260.

- J1. Ambulance battery 12v CC
- J2. Battery connection
- J3. Light signal LED right
- J4. Light signal LED left
- J5. Electro magnet (Brake)
- J6. Inverter Motor connection line (24v CC)
- J7. Micro Switch
- J8. Motor hall signal
- J9. Motor control
- J10. Terminal connection Control Panel card to Tocuh Panel (D2)


06 OPERATION

REV. 2026/07

ELECTRICAL WIRING DIAGRAM:

The below diagram shows the electrical scheme of the track PA-260.

- A. Power input line (12v CC)
- B. Microswitch line load terminal
- C. PCB control panel line to main PCB
- D. Motor 200w connection line (24v CC)
- E. Touch panel connection to PCB control panel
- F. Supply line left LED lights
- G. Supply line right LED lights
- H. Load line between support and ambulance (12v CC)
- I. Battery connection to main PCB

06 OPERATION

REV. 2026/07

PROBLEM SOLVING:

COD.	DESCRIPTION OF THE ERROR	POSSIBLE CAUSE	POSSIBLE SOLUTIONS
1	- THE BATTERY VOLTAGE IS LESS THAN 3V - IS DISCONNECTED OR DAMAGED	- BATTERY DISCONNECTED - BATTERY DAMAGED	- CONNECTING THE BATTERY CORRECTLY - CHECK CONNECTION OF THE BATTERY CABLES TO THE ELECTRONIC PLATE - REPLACE BATTERY
2	- THE BATTERY VOLTAGE IS SUPERIOR TO 34V	- OVERLOAD DURING BATTERY CHARGING	- COOL THE BATTERY FOR 3/4 HOURS - USING THE CHARGER PROVIDED BY PROMEBA - REPLACE BATTERY FOR A NEW ONE
3	- THE STRETCHER DOES NOT RECEIVE TENSION FROM AMBULANCE	- FAILURE TO CONNECT SUPPORT-TRACK - AMBULANCE BATTERY DOWNLOADED OR DAMAGED - POWER SUPPLY CABLES CONNECTED BETWEEN THE SUPPORT AND THE AMBULANCE - ONE OF THE COPPER TERMINALS OF THE TRACK HAS NOT CORRECTLY LEFT	- CHECK CONNECTION TERMINALS BETWEEN THE CATERPILLAR AND THE SUPPORT - CHARGING BATTERY OR REPLACING FOR A NEW ONE - CHECK CONNECTION OF THE CABLES BETWEEN ORUGA-SUPPORT-AMBULANCE - CHECK THAT THE CONNECTION TERMINALS OF THE TRACK ENTER AND EXIT FROM YOUR ACCOMMODATION CORRECTLY
4	- THE BATTERY VOLTAGE IS LESS THAN 18V	- THE BATTERY IS DOWNLOADED	- LOADING BATTERY TO THE MAXIMUM LEVEL
5	- THE TRACK RECEIVES TENSION THROUGH THE TERMINAL CONNECTION TERMINALS SUPPORT BUT NOT HAVE BEEN PRESSED COMPLETELY	- TRACK HAS NOT ENTERED PROPERLY - MICROPOWER FAILURE	- MAKE SURE THAT THE TRACK HAS BEEN FITTED INTO THE SUPPORT AND THAT THE CONNECTING TERMINALS ARE WELL CONNECTED - CHECK MICROPOWER STATUS
6	- BATTERY IS NOT DETECTED	- BATTERY DISCONNECTED - CONNECTING THE DAMAGED BATTERY - VERY LOW BATTERY LOAD	- ADJUST THE BATTERY TO YOUR BASE - REPLACE FOR NEW BATTERY - CHARGING BATTERY WITH EXTERNAL CHARGER
7	- FAILURE OF ENGINE ENCODER	- CONTACT WITH TECHNICAL SERVICE	- CONTACT WITH TECHNICAL SERVICE
8	- INTENSIDAD SUPERIOR A 48V	- TRACK WITH WEIGHT OVERLOAD - LOW BATTERY - THE TRACK DOES NOT UP OR LOWER SOFTLY	- DO NOT LOAD TRACK WITH MORE THAN 180KG - CHARGING BATTERY WITH EXTERNAL CHARGER - LOOK FOR POSSIBLE HITS, INTERFERENCES OR PARTS LOCKED IN THE TRACK
9	- INTENSITY UPPER 48V	- TRACK WITH WEIGHT OVERLOAD - LOW BATTERY - THE TRACK DOES NOT UP OR DOWN WITH SMOOTHNESS	- DO NOT LOAD TRACK WITH MORE THAN 180KG - CHARGING BATTERY WITH EXTERNAL CHARGER - LOOK FOR POSSIBLE HITS, INTERFERENCES OR PARTS LOCKED IN THE TRACK
10	- FAILURE OF THE LOADING PROCESS	- CONTACT WITH TECHNICAL SERVICE	- CONTACT WITH TECHNICAL SERVICE

IMPORTANT CONSIDERATIONS

The electronic track system will automatically shut down if no operation has been performed in 15 minutes.

In up and down maneuvers repeatedly of up to 10 cycles uninterruptedly a 15 seconds stop must be made mandatory for recovering the motor-reducer system of the track.

It is advisable to remove the battery from the interior of the track if it is going to go through long periods of non-activity.

It is important to maintain a minimum of charge cycles every month to maintain good battery capacity. Do not leave the battery uncharged for more than 2 MONTHS. It could stop working.

07 GENERAL MAINTENANCE

REV. 2026/07

07.1. Cleaning

It is essential to keep the equipment clean to ensure proper use and durability of the set. Thorough cleaning should be carried out periodically, especially in areas exposed to dirt that can accumulate dirt such as corners or holes.

Do not use high pressure cleaning systems or bleach to disinfect the product, they can damage it. Instead, use a hydroalcoholic-based disinfectant, wash it with water and let it dry naturally, do not use direct heat sources to dry.

Non-compliance with cleaning operations can lead to the risk of cross infections due to the presence of secretions and/or residues. During all control and disinfection operations, the operator must wear appropriate personal protective equipment, such as gloves, glasses, etc.

When the cleaning procedure includes a disinfectant product, carefully follow the instructions of the manufacturer of the product used regarding the method of application and contact time. Make sure that you have taken all the necessary precautions to ensure that there are no risks of cross-infection or contamination of patients and operators.

Failure to carry out the correct cleaning routine could increase the risk of cross infection, due to presence of body fluids and/or residuals.
The operator must always wear adequate personal protection such as gloves and mask, etc. during all checking and cleaning procedures.

Be sure to keep these zones free of dirt.
Have special care in the cleaning to ensure good functioning of the set

07.2. Lubricating

In general, all moving parts must be greased. Our products leave the factory completely greased and lubricated, however, it is possible that over time and use of the product the elements become degreased, either due to loss of lubrication or due to dirt.

Periodically clean and oil affected areas according to manufacturer's specifications. Check for loose, missing or worn parts.

Never lubricate the belts. Lubrication can cause the belts to perform unpredictably, which can cause injury to the patient and/or operators.

07.3. Wear areas

Periodic inspection of the different components of the system for signs of wear is a preventative measure that can reduce breakdowns. Check possible lubricant leakages,

07 GENERAL MAINTENANCE

REV. 2026/07

grooves or bearing in poor condition. Check also, due to intensive usage, the moving parts such the handles, the wheels and the folding system of the chair and the track.

In the event that a worn component needs to be changed, you should contact our sales department for more information on ordering spare parts and their installation.

The repair or replacement of product components must necessarily be carried out by qualified personnel, that is, who have received training from the Manufacturer, Promeba S.L. and original spare parts should always be used.

Promeba S.L. declines any responsibility for any damage, direct or indirect, resulting from improper use of spare parts and/or any repair intervention carried out by unauthorized persons.

07.4. Maintenance program

Although the frequency of the controls is determined by factors such as legal requirements, the type of use, the frequency of use, the environmental conditions during use and storage, it is placed a recommendation on the maintenance that the equipment must have for correct operation throughout its useful life.

All maintenance and revision activities must be recorded and documented with the corresponding technical intervention reports, and the documentation must be kept for at least 10 years from the end of the product's useful life and must be made available to the competent Authorities and/or from the Manufacturer, when requested.

Maintenance review	Every use	When necessary	Every month	Every year	Every 4 years
Disinfect	X				
Clean		X			
Inspect	X	X			
Grease necessary components		X		X	
Change springs				X	
Change wheels					X

For both the disinfection and cleaning procedures, point 7.1 of this manual must be followed.

During the inspection procedure, which must be done after each use, the following tasks must be carried out:

- Check that all components are present.
- Check the integrity of the device, that is, that there are no breaks, cracks, holes or cuts.
- Check the wear of the product following section 7.3 of this user guide.
- Check that the moving parts slide correctly.
- Check that the accessories provided are present, intact and that they function correctly.

For some conditions of use, due to vibrations or impacts, certain elements may lose their tightening torque or fixing properties. Periodically check for loose items, especially moving parts. Inspect any mechanical connections. Always observe and adhere to the recommended tightening torques.

08 SPARE PARTS


REV. 2026/07

To request a spare part, indicate:

- Number and description of the set in which the necessary part is found (see table)
- Point out the piece in the drawing, or the area where it is located.

Following these instructions, our sales department will provide you with a detailed drawing of the assembly so that you can specify the reference of the required spare part.

PS2510-00000

08 SPARE PARTS


REV. 2026/07

Nº	COMPONENTS DESCRIPTION	UNITS
1	PS2500-00041 - Casquell gir potes posteriors cadira	2
2	PS2500-00051 - CASQUELL GIR MANETEAS RESPATLLER CADIRA	2
3	PS2500-00110 - SEIENT CADIRA	1
4	PS2500-00130 - Conjunt cinturó cadira PS-250	1
5	PS2500-00130 - CONJUNT CINTURÓ RESPATLLER CADIRA PS-250	1
6	PS2510-01000 - CONJUNT ESTRUCTURA RESPATLLER CADIRA	1
7	PS2510-02000 - CONJUNT POTA DAVANTERA CADIRA	1
8	PS2510-04000 - CONJUNTO MANETA EXTENSIBLE CADIRA	2
9	PS2510-04071 - CASQUELL GIR MANETES EXTENSIBLES CADIRA	2
10	PS2510-06002 - SUB. SOLDADURA REPOSAPÉUS CADIRA	1
11	102032059901 - ARANDELA CURVA TUB 30-34MM	4
12	PS2500-00122 - RESPATLLER CADIRA	1

08 SPARE PARTS

REV. 2026/07

PA2600-00000

08 SPARE PARTS



REV. 2026/07

Nº	COMPONENTS DESCRIPTION	UNITS
1	PA2600-01002 - SUBCONJUNT ESTRUCTURA MARC ORUGA	1
2	PA2600-02002 - SUBCONJUNT ESTRUCTURA ORUGA	1
3	PA2600-03000 - CONJUNT QUADRE DE MANDOS EXTENSIBLE	1
4	PA2600-04002 - SUBCONJUNT REGULACIÓ DESPLEGAT ORUGA	1
5	GN 822-7-B-ST - POM POSICIONADOR SENSE BLOQUEIG	1
6	PA2600-07000 - Conjunt LED encastat oruga	2
7	PA2600-05000 - CONJUNT CAPÇAL DE CARREGA ORUGA	1
8	696Z - Rodament de boles sellat 5x15x6mm	4
9	626ZZ - RODAMENT DE BOLES SELLAT 6X19X6MM	4
10	PA2600-00011 - CASQUELL TENSOR CINTA ORUGA	2
11	PA2600-00030 - Politja dentada oruga D56	2
12	PA2600-00040 - POLITJA DENTADA ORUGA D85	4
13	PA2600-00050 - Eix gir conjunt oruga	1
14	F6801Z - Rodament de boles sellat 12x21x5mm	2
15	PA2600-00060 - Plat unió politges	1
16	PA2600-00080 - Corretga dentada oruga	2
17	TAP-PLASTØ30 - TAPON PLASTICO FINAL Ø30 - Mod. 17-30	1
18	PA2600-00101 - Tub regulació interior desplegat oruga	1
19	PS2500-00031 - BALDA BLOQUEIG PLEGAT CADIRA	2
20	PC6500-03021 - molla compressió	2
21	PA2600-00140 - CASQUELL REGULADOR INCLINACIÓ ORUGA	1
22	PA2600-00172 - TAPA MANETAS BLOQUEIG MARC EXTENSIBLE	1
23	PA2600-00192 - Tivador desbloqueig capçal extensible	2
24	PC7500-02100 - BATERIA ENERGIVM 2500 MAH - 25.6V	1
25	PC7500-02110 - Adaptador bateria EnergiVm	1
26	PA2600-00090 - Plat dentat motor oruga	1
27	PA2600-00130 - Pinyo dentat motor oruga	1
28	PA2600-00072 - TENSOR CORRETGA ORUGA ESQUERRA	2
28	PA2600-00022 - Tensor corretga oruga dreta	2
30	PA2600-00110 - PLACA DE CONTROL PRINCIPAL ORUGA PA-2601	1
31	CRLCBSRE-4-01BK - SEPARADOR PCB H6.4	4
32	PA2600-00020 - Motor KELVIN K80-63125	1
33	PA2600-00202 - Safata suport bateria oruga	1
34	PA2600-00210 - Conector DC encastable 2.1mm femella	1
35	PB4400-00410 - Pulsador Rocker (RL23-23-C-6-BKBK-P4-3)	1
36	PA2600-00242 - Xapa desloqueig balda oruga	2
37	PA2600-00252 - Tapa manetas bloqueig oruga	1
38	PA2600-01020 - Tub ganxo fixació marc oruga	2
39	3069 - RUEDA DIAM.50MM REF. ARSA - 3069	2

08 SPARE PARTS

REV. 2026/07

PA2810-00000

Nº	COMPONENTS DESCRIPTION	UNITS
1	A-8-12-20 - CASQUELL BRONZE 8X12X20	1
2	PA2800-00012 - XAPA GANXO SUPERIOR	1
3	PA2800-00020 - PROTECTOR PLASTIC INFERIOR LATERALS	2
4	PA2800-00030 - PROTECTOR PLASTIC SUPERIOR LATERALS	2
5	PA2810-01002 - Sub. soldadura estructura cos inferior	1
6	PA5700-00041 - Molla PA570	1
7	TAP-PLASTØ25 - TAPON PLASTICO FINAL Ø25 - MOD. 17-25	2

REV. 2026/07

Keep this document for at least 10 years after the end of the device's useful life.

102

REV. 2026/07

Keep this document for at least 10 years after the end of the device's useful life.

103

11 LEGAL NOTICES

REV. 2026/07

- This document may contain technical inaccuracies or typographical mistakes.
- Changes are periodically added to the information herein; these changes will be incorporated in new editions of the publication.
- Promeba, S.L. reserves the right to make any modification or improvement in the products described in this publication if it is appropriate.
- Promeba, S.L. may have patents or patent applications which address themes described in this document. The possession of this document does not entitle no license to these patents. The information contained in this document does not affect or change the specifications or product specifications or warranties of Promeba, S.L.
- No part of this document shall operate as express or implied license or compensation under the intellectual property rights of Promeba, S.L. or .third party.
- All the information contained in this document has been obtained in specific environments and it is presented as an illustration. The results obtained in other operating environments may vary.
- Promeba, S.L. can use or distribute the information that the customer provides in whatever way they see fit, without incurring any obligation with the client.

12 PRODUCT WARRANTY

REV. 2026/07

Promeba, S.L. guarantees that its products have satisfactorily passed all established quality controls, both functional and material. The duration of the guarantee is 2 years from the date of purchase of the product.

This guarantee will be valid only when it is presented with the original invoice or proof of purchase (indicating the date of purchase, model and the name of the distributor) together with the defective product during the period covered by the guarantee. Promeba, S.L. reserves the right not to offer the free warranty service if the indicated documents are not presented or if the information they contain is incomplete or illegible.

1. This warranty will not apply if the model name or serial number of the product has been altered, erased, disappeared or becomes illegible.

2. This guarantee does not cover transportation costs or risks derived from the transportation of your product to and from Promeba, S.L.

3. This warranty does not cover any of the following:

- a) Periodic maintenance and repair or replacement of parts derived from normal wear and tear.
- b) Expendable material (components that are expected to need periodic replacement during the life of the product, such as belts, straps, etc.).
- c) Damages or defects derived from the improper use, operation or treatment of the product and not due to normal use of the product.
- d) Damages derived from:
 - i. Misuse, including:
 - Treatment that results in damages or physical, superficial or

appearance changes of the product.

- Installation, use or storage of the product in a way that does not respect the instructions described by Promeba, S.L.

- Maintenance of the product in a way that does not respect the instructions of Promeba, S.L. for proper maintenance.

- Installation or use of the product in a way that does not respect the technical or safety regulations of the country where it is used or installed.

ii. Use of components not provided with the product or incorrect installation of accessory parts not previously tested.

iii. States or defects of the system in which the product is used or incorporated with the exception of other products Promeba, S.L. designed for use with the product.

iv. Use of the product with accessories, peripheral units and other products of a type, condition or standards not established by Promeba, S.L.

v. The manufacturer or distributor will be solely responsible for deciding whether to send the parts for repair, or to replace the product in its entirety. In no case will operators be sent to repair or replacement of the product.

Except in the cases mentioned above, Promeba, S.L. will make no warranties in relation to product, performance, accuracy, reliability, or fitness for logic purpose of the equipment or other type. If this exception is not legal or contemplated by current law, Promeba, S.L. will limit or exclude your warranties only to the extent permitted by applicable law.

12 PRODUCT WARRANTY

REV. 2026/07

The only obligation on the part of Promeba, S.L. in connection with this warranty is to repair or replace parts subject to the terms and conditions of this warranty.

Promeba, S.L. is not responsible for loss or damage to products, this warranty or others, including economic loss or non-assessable damage; the price paid for the product; loss of profits, income, information, usufruct or use of the product or associated products or indirect, accidental or critical loss or damage.

This clause refers to whether the loss or damage is due to deterioration or inoperability of the associated product due to defects or unavailability of Promeba,

S.L., which has caused downtime, loss of user time or business interruption.

In cases where the law prohibits or limits these liability exclusions, Promeba, S.L. will exclude or limit his liability only to the extent permitted by applicable law.

For example, there are countries that prohibit the exclusion or limitation of damages caused by negligence, reckless negligence, willful misconduct, fraud and similar actions. The responsibility of Promeba, S.L. In this guarantee will not exceed, in any case, the price paid for the product, but if the current law allows only limitations of greater responsibilities, these will apply.

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*All rights reserved. Variations can be done whitout notice.