MANUAL DE INSTRUCCIONES // USER GUIDE
ELEVADOR CORTO PB-440
ELECTRIC STRETCHER HOLDER PB-440 SHORT

DT14_02_05_REV. 2026/03

00 INDEX

REV. 2026/03

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00 INDEX

REV. 2026/03

01	INTRODUCTION.....	36
01.1	Products.....	36
01.1	Using this guide.....	36
01.2	Legend of symbols.....	36
01.3	Service request.....	36
01.4	Life span.....	37
01.5	Demolition.....	37
01.6	Reference standards.....	38
01.7	Labelling.....	39
02	INTENDED USE.....	40
02.1	Intended use.....	40
02.2	Contraindications and side effects.....	40
02.3	Physical requirements of the operators.....	40
03	WARNINGS.....	41
03.1	General warnings.....	41
03.2	Specific warnings.....	43
03.3	Residual risk.....	44
04	DESCRIPTION OF THE PRODUCT.....	45
04.1	Main components.....	45
05.2	Technical data sheet.....	46
05	OPERATION.....	47
05.1	Transport and storage.....	47
05.2	Electric connection.....	47
05.3	Anchoring points.....	49
05.4	Operation.....	49
05.5	Electric movements inside the ambulance.....	53
05.6	Manual operation of the support.....	54
05.7	Electrical diagram.....	55
06	GENERAL MAINTENANCE.....	58
06.1	Cleaning.....	58
06.2	Lubrication.....	58
06.3	Wear areas.....	59
06.4	Maintenance program.....	59
07	SPARE PARTS.....	60
08	TRAINING REGISTER.....	61
09	MAINTENANCE REGISTER.....	62
10	LEGAL NOTICES.....	63
11	PRODUCT WARRANTY.....	64

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REV. 2026/03

01.1. Products

REFERENCE	UDI CODE	TRADE NAME
PB4402-00000W	8436601872067	ELECTRIC STRETCHER HOLDER PB-440 SHORT

01.2. Using this guide

The manual provides instructions for use and maintenance of the product, as well as technical aspects, operation, spare parts and safety.


It is recommended before the operation of the product to read carefully this manual in order to avoid damages caused by a misuse, regardless of the level of experience acquired in the past with similar devices.

Do not lose this document. It should be accessible to any doubt that could appear by medical personnel. In case the device has some instructions for use of a different product than the one received, it is necessary to contact the manufacturer immediately before using the product.

Remember that a good use and maintenance are necessary for the proper operation of the product.



Each product includes an identifying adhesive label with the serial number and model. Keep these numbers so that you can refer them to the dealer if necessary.

01.3. Legend of symbols

SYMBOL	EXPLANATION / DESCRIPTION
	MANUFACTURER symbol. This symbol is accompanied by the name and address of the manufacturer, adjacent to the symbol (PRODUCTOS METÁLICOS DEL BAGÉS S.L., Ctra. C-16 Km 59.5, 08650 Sallent (Barcelona)).
	Indicates the manufacturer's reference number to identify the medical device. PROMEBA, S.L. uses this symbol to set each internal reference for each configuration and business variant
	Indicates the manufacturer's serial number to identify a specific medical device.
	Indicates the manufacturing date. The symbol must be accompanied by a manufacturing date (yyyy-mm), adjacent to the symbol.
	It is placed to inform that the product is a "Medical Device"
	CE symbol without the intervention of a Notified Organism, as a Medical Device classified as Class I according to the EU Regulation 2017/745 on Medical Devices.
	Symbol for the Unique Device Identifier.
	Symbol "See instructions for use or operating instructions"

01 INTRODUCTION

REV. 2026/03

	Symbol "caution". This symbol is placed to warn of the need for the user to refer to important precautionary information in the operating instructions, such as warnings and cautions not otherwise found on the label.
	Waste symbol for electrical and electronic equipment "WEEE" (Waste Electrical and Electronic Equipment). Recycle: Electronic equipment. DO NOT THROW GARBAGE
	"Dangerous voltage" symbol. To warn of electricity.
	Symbol "Caution". For a general warning.
	Warning, crushing of hands
	"Double insulation" symbol. The protection of appliances marked with this symbol is ensured by double insulation and does not require a safety electrical connection to earth (ground).
	Cot duty cycle: 16.7% (less than 60 seconds on, more than 300 seconds off)
	Indicates the temperature limits. The upper and lower temperature limits should be indicated adjacent to the horizontal lines.

01.4. Service request

For information of the correct interpretation of the instruction manual, the use, maintenance, installation and restoration of the product, please contact Promeba customer service: T. 93 837 12 00, email: promeba@promeba.com or write to PROMEBA, S.L. - Ctra C-16 Km 59.5 - 08650 Sallent (Barcelona) - SPAIN.

Whenever you want to communicate any matter related to a product, you will be asked for the product reference (**REF**) and the serial number (**SN**). It is strongly recommended to keep the label containing this information and found on the original packaging of the product.

01.5. Life span

If used as described in the following instructions, this device has a useful life of 10 years from the date of purchase.

This useful life can be extended with annual reviews carried out by the manufacturer, which uses specialized and authorized internal and external technicians.

In case these annual checks are not carried out, the device must be disposed of according to the information in paragraph 02.5 and the manufacturer must be notified.

Only the manufacturer or an authorized center can extend the life of the device, if it meets the safety and legal requirements.

Promeba, S.L. will not accept any responsibility for malfunction or damage caused by the use of devices that have not been checked by the manufacturer or authorized center, or that have exceeded the maximum allowed useful life.

01 INTRODUCTION

REV. 2026/03

01.6. Demolition

When the devices are no longer suitable for use, they must be correctly recycled according to whether it is electrical and electronic equipment, batteries or normal solid waste, following current regulations on demolition:

WARNINGS FOR THE CORRECT DISPOSAL OF THE PRODUCT IN ACCORDANCE WITH EC DIRECTIVE 2012/19/EU WEEE:

At the end of his life, the product must not be disposed as household waste. Can be taken to special recycling centers provided by local government, or return it to the dealer on purchase of a new device of the same type and used for the same functions. Dispose of the product separately avoids possible negative consequences for the environment and human health resulting from inappropriate disposal and allows to recover the materials in order to obtain significant savings in energy and resources. The symbol on the label indicates separate collection of electrical and electronic equipment.

An incorrect disposal of electrical and electronic equipment could result in sanctions.



01.7. Reference standards


REFERENCE	TITLE OF DOCUMENT
EU Regulation 2017/745	EU Regulation on Medical Devices
EN 1789:2007+A2:2014	Medical vehicles and their equipment - Road ambulances
UNE-EN 1865-5:2012	Patient handling equipment used in road ambulances - Part 5: Stretcher support
EN 60601-1:2006/A1:2013 IEC 60601-1:2005/A1:2012	Medical electrical equipment - Part 1: General requirements for basic safety and essential performance
EN 61000-6-1:2007 IEC 61000-6-1:2005	Electromagnetic compatibility (EMC) - Part 6-1: Generic standards - Immunity for residential, commercial and light-industrial environments
EN 61000-6-3:2007+A1:201 2 IEC 61000-6-3:2006+A1:201 0	Electromagnetic compatibility (EMC) -- Part 6-3: Generic standards - Emission standard for residential, commercial and light-industrial environments

As a distributor or end user of the products manufactured and/or distributed by Promeba, S.L., you are strictly obliged to know the legal provisions in force in the country of destination of the goods, applicable to the devices to be supplied (including the regulations relating to technical specifications and/or safety requirements) and, therefore, understand the necessary requirements to ensure compliance of the products themselves with all the legal requirements of the territory.

01.8. Labelling

Each product incorporates an identification label, placed on the device itself. It must never be removed or covered. This label includes the data of the company which produced the device (), the product code / reference (), the manufacturing date (), the serial number (), the UDI-PI code () and the CE distinctive (). Please keep these numbers so you can inform the dealer if necessary.

The UDI-PI code is made up of the following fields:

- (01)08436601872067: UDI code of the product.
- (11)260321: it gives information about the batch release date of the product following the format *aammdd*, where *aa* is the batch release year, *mm* is the month of the batch release and *dd* is the day of the batch release.
- (21)SN2600000: identifies the serial number of the product giving information about the batch release year (26) and a correlative number following the manufacturing order *nnnnn* (00000).

Each product incorporates an identification label, placed on the package.

REF PB4402-00000W	
 2026-03	
	SN2600000 SN
ELEVADOR PB-440 CORTO ELECTRIC STRETCHER HOLDER PB440 SHORT PORTE-BRANCARD ÉLECTRIQUE PB440 COURT ELEKTRISCHER TRAGSTELLER-HALTER PB440 KURZ PORTA BARELLA ELETTRICA PB440 CORTO	
  	
	
UDI	(01)08436601872067(11)260321(21)SN2600000
 Productos Metálicos del Bages, S.L. Ctra. C-16 km. 59,5 08650 - SALLENT BARCELONA (SPAIN)	MADE IN E.U. FOR-181.Rev. 2025/09

02 INTENDED USE

REV. 2026/03

02.1. Intended use

The product ELECTRIC STRETCHER HOLDER PB440 SHORT is a medical device intended for installation in ambulances to support, secure, and electrically assist the loading, unloading, elevation, and positioning of stretchers carrying patients during medical transport, minimizing manual handling by healthcare personnel.

02.2. Contraindications and side effects

Any contraindication or side effect have been resulted from the clinical analysis of the ELECTRIC STRETCHER HOLDER PB440 SHORT family of products.

02.3. Physical requirements of the operators

Promeba stretcher supports are intended for professional use only. The device is intended to be used by medical transport technicians or ambulance technicians who have knowledge and experience with patient transport.

Operators must be able to assess the integrity of the fixing, knowing the mechanism for fixing the stretcher to the support and being able to ensure its correct operation. They must also be able to assess any anomaly in the system/assembly by communicating the problem to those responsible and stopping the use of the devices if necessary to ensure patient safety.

The capacities of the various operators must be considered before determining his role in the employment of the board.

03 WARNINGS

REV. 2026/03

03.1. General warnings

1. The product must be used by trained personnel only, having attended specific training for this device and not for similar products.
2. Establish a maintenance and periodic control program and identify a reference manager. The person in charge of the ordinary maintenance of the device must guarantee the basic requirements provided by the manufacturer in these instructions for use. It is recommended that the subject check at least every 6 months if there are updated user manuals available or if there are any changes related to his / her product. This information is freely available on the website <http://promeba.com/>.
3. Training routines must be registered on a special register in which the names of those trained, of the trainers, date and place are indicated. This register, which will certify the eligibility of the operators to use the Promeba, S.L. device, has to be kept for a period of 10 years after the disposal of the device itself. This register will be made available to the competent authorities and/or manufacturer if requested.
4. Promeba, S.L. is always at your disposal to plan trainings on products.
5. Before carrying out any kind of operation on the appliance (training, installation, use), the operator must carefully read the enclosed instructions, paying particular attention to the correct safety precautions and to the procedures to be followed for installation and for correct use.
6. If the instructions belong to another device and not the device received, inform the manufacturer immediately and avoid use of the device.
7. In the case of any doubts as to the correct interpretation of the instructions, please contact Promeba, S.L. for any necessary clarifications.
8. Do not allow untrained persons to help during the use of the device, because they could cause damage to the patient or to themselves.
9. Before each use of device the perfect operating state of the device must be checked as specified in the instruction manual. If any damage or abnormalities which could in any way influence the correct functioning and the safety of the device, of the patient and/or of the user are detected, the device must be immediately removed from service and the manufacturer must be contacted.
10. If any failure or incorrect functioning of the device is detected, it must be immediately substituted with a similar item so that the rescue procedures are guaranteed without any interruption.
11. Use of the device in anyway other than described in this manual is forbidden.
12. Do not alter or modify in any way the device; any such interference could cause malfunctions and injury to the patient and/or rescuer.
13. The device must not in any way be tampered with (modification, adjustment, addition, replacement). In such cases all responsibility will be denied for any malfunctions or injuries caused by the device itself; moreover CE certification and product warranty will be considered void.

03 WARNINGS

REV. 2026/03

14. Use only original or PROMEBA approved components/spare parts and/or accessories to carry out all operations without causing alterations or modifications to the device; Otherwise, any responsibility for the operation is declined.
15. Ensure that all the necessary precautions are taken in order to avoid the hazards that can arise as the result of contact with blood or body fluids.
16. Register and store with these instructions: lot number, place and date of purchase, first date of use, date of checks, name of users and any other comments.
17. When the device is being used, the assistance of qualified staff must be guaranteed.
18. The product must not be exposed or come into contact with thermal sources of combustion and flammable agents, but must be stored in a dry, cool place, protected from light and sunlight.
19. Store and transport device in its original packaging.
20. Do not store the product under other more or less heavy materials, which can damage the structure of the product.
21. Attention: laboratory testing, post production tests and instruction manuals cannot always consider every possible scenario for use. This means that in some cases the performance of the product could be notable different from results to date obtained. Instructions are continually being updated and are under tight surveillance of fully qualified staffs with adequate technical formation.
22. As a distributor or end users of products manufactured and/or marketed by Promeba, S.L., you are strictly required to have a basic knowledge of any legal requirements existing in the goods final destination country applying to the devices contained in this supply (including laws and norms regarding technical specifications and/or safety requirements) and therefore you are also strictly required to have the necessary knowledge to guarantee all aspects regarding the total conformity of the products to the regulations in the relevant territory.
23. Promptly notify PROMEBA regarding any revisions to be made by manufacturer in order to guarantee the conformity of the product to the territory's legal specifications (including those resulting from rules and/or norms of other nature).
24. Act with all due care and diligence, and contribute to ensure conformity to general safety requirements of all devices marketed in the territory, by providing final users with all necessary information for carrying out periodical checks on their devices, as specified in the relevant user manual.
25. Actively contribute to product safety checks on products sold, by communicating any relevant risk analysis information both to the manufacturer and to any competent authorities so that the necessary actions can be promptly taken.
26. Be aware that in the event of any failure to conform to the above mentioned requirements you will be deemed fully responsible for all damages that might occur. Therefore we expressly disclaim any responsibility and/or liability for your non-compliance with the present "Regulatory provisions".
27. Follow internal procedures and protocols approved by your organization.

03 WARNINGS

REV. 2026/03

28. Follow approved Emergency Medical Service procedures for patient immobilization and transportation.

03.2. Specific warnings

1. Use only accessories/spare parts that are original or approved by Promeba, S.L. when performing any operation so as not to cause any alteration or modification to the device, otherwise we do not assume any responsibility for the correct operation or resulting damage of the device to the patient or the operator and the guarantee will be considered void, as well as compliance with the Regulation 2017/745.
2. Always respect the maximum capacity of the device, as indicated in this user's manual. Maximum load capacity means the total weight distributed according to the human anatomy. In determining the load of the total weight on the product, the operator must consider the weight of the patient, the equipment and the accessories. Moreover, the operator must consider that the overall dimensions of the patient do not reduce the functionality of the device.
3. Inside the ambulance, place and adjust the device taking care not to obstruct other parts or other equipment. Also consider lateral displacement if it is available in your unit.
4. Never leave the patient unassisted on the device, because he may be injured.
5. Do not use bleach to disinfect the device. Use a hydroalcoholic based disinfectant and wash with water.
6. The device and all its components, after washing, should be allowed to dry completely before storing.
7. Avoid contact with sharp objects.
8. Do not use the device if it is pierced, deformed or shows signs of loss of integrity.
9. Avoid dragging rough objects on the surface of the device.
10. Make sure the spine board or scoop stretcher is securely secured by using the rail belts before moving the ambulance. If it is not properly fixed, there is a risk of impact.
11. Do not operate in case the weight has not been distributed correctly.
12. Condensation, water, ice and dust accumulations can affect the correct operation of the device, making it unpredictable.
13. Always immobilize the patient at least with the belts supplied by the manufacturer during any movement of the stretcher, including loading and unloading to the ambulance.
14. Do not use the device as a support for climbing or leaning on and reaching high points in the ambulance. The product may be damaged and the technician or patient may be harmed.
15. Disinfection activities must be carried out in accordance with the validated cycle parameters, indicated in the specific technical standards.

03 WARNINGS

REV. 2026/03

16. The stretcher used together with the stretcher support, must be installed in accordance with the specifications provided by the manufacturer and ensuring its compatibility with technical specifications of the stretcher support.
17. Promeba S.L. shall not be liable if unauthorized / incompatible stretcher fastener is used.
18. The installation of the stretcher support on the vehicle floor must be carried out ensuring compliance with the safety requirements of EN 1789. In order to be in compliance with the requirements of EN 1789, the installer has the duty to check the fixation based on the exact position on the floor and based on the structure of the floor itself.
19. Do not place body parts under the stretcher support; this may cause crushing risks when the stretcher support is moved. Never put your feet or hands under the device.
20. Installation carried out by untrained personnel, could result in inadequate fastening of the device inside the vehicle, with consequent risks related to its instability or mobility.

03.3. Residual risk

The residual risks listed below have been identified only with reference to the intended use of the device:

1. Do not place body parts under the stretcher support; this may cause crushing risks when the stretcher support is moved. Never put your feet or hands under the device.
2. Installation carried out by untrained personnel, could result in inadequate fastening of the device inside the vehicle, with consequent risks related to its instability or mobility.
3. Use by untrained personnel may result in injury to patient, operator, or third parties.


04 PRODUCT DESCRIPTION

REV. 2026/03

04.1. Main components

PB4402-00000 ELECTRIC STRETCHER HOLDER PB440 SHORT

N°	DESCRIPTION OF THE COMPONENTS	N°	DESCRIPTION OF THE COMPONENTS
1	Crank	7	Emergency stop pushbutton
2	Third anchor point	8	Side movement pedal
3	Security sensors	9	Trendelemburg position button
4	Extraction lever	10	Status indicator led
5	Manual loading lever	11	Manual unlocking hooks
6	Charge / discharge button	12	On-off button





04 PRODUCT DESCRIPTION

REV. 2026/03

04.2. Technical specifications / data sheet

Models	PB4402-00000W
Length	2018mm
Width	630mm
Height	425mm
Weight	112kg
Max. Load	300kg
Electric supply	12V CC
Power consumption	100 ^a
Trendelemburg	-15° / +10°
Height loading	± 315 mm
Side Movement	300mm

05 OPERATION


REV. 2026/03

05.1. Transport and storage

Before transporting the device, make sure that it is correctly packaged.

Transport with a crane (high load)


Thread tapes through the upper platform for its transport (as shown in graphs). Transport the load levelled following all the precepts and regulations for the transport of suspended loads.

TRANSPORT WITH A CRANE
FIX THE TAPES TO LEVEL THE LOAD WEIGHT

Transport with forklift

Place a pallet or flat surface under the rail and transport transversally. Do not load the rail directly on the forklift without a pallet or flat support surface. Do not load the rail longitudinally on the forklift if you do not have a sufficiently long base. Each rail has been thoroughly inspected upon leaving the factory. To ensure that it has not been damaged during transport, please carefully inspect the exterior and interior, and if you find any damage, immediately notify the Promeba or the authorized distributor. For a correct transport operation, the rail box must be levelled.

TRANSPORT WITH A FORKLIFT
USE ALWAYS A FLAT SUPPORT SURFACE

To unpack the device, place the box on a flat, stable surface and carefully remove the packing tape. Remove the device from the inside of the case carefully to prevent shocks to the device. Please retain the original packaging for use in case of any other transportation and for storage.

Product damage caused during shipping and handling is not covered by the warranty. Repair or replacement of damaged parts is the responsibility of the customer.

The device should be stored in a cool, dry place, away from direct sunlight. It must not be in contact with substances or chemical agents that may cause damage and reduce its safety performance.

The storage temperature must be between -20°C to 40°C.

The product must not be manipulated in any way (modification, adjustment, addition, substitution). In such cases, all responsibility for any malfunction or injury caused by the product itself will be disclaimed; furthermore, CE certification and product warranty will be considered void.


05.2. Electric connection

To make the initial connection, follow the steps described below.

Be especially careful to connect the electric stretcher holder to the required 12V DC power supply and do not leave the cables in high-traffic areas or near moving parts, as this could cause breaks, cuts, and system failures.

05 OPERATION

REV. 2026/03

B: To make the connection, look for the black cable placed in the front of the equipment. Strip the plastic insulation from the ends of the red and black wires using a wire stripper so you can make the connection.

C: Use a power strip and a flathead screwdriver to connect the red and black wires to those on the ambulance's power outlet.

MAKE SURE THE POWER LINE IS PROPERLY PROTECTED AND THAT THE POWER SUPPLY ALONG THE LINE IS STABLE.


05 OPERATION

REV. 2026/03

05.3. Anchoring points

Before putting the electric stretcher holder PB-440 short into operation, you must ensure that each and every anchor point is in place and securely fastened.

There are 8 anchoring points, as described in the image below.

05.4. Operation

LOADING AND UNLOADING OPERATION

Starting from the position in which the bench has been left after the unloading (A):

LOADING

- Push the ON button on the rear part of the electric support (A).
- Remove the upper part of the support by pressing the “EXTRACTION” handle (see Figure 1) and pulling it out of the ambulance (see Figure A).
- Using the button located on the right side of the upper part of the support (see image 2), press (I) to perform the AUTOMATIC LOADING operation; wait until the upper part of the support is fully tilted (see image B).
- The locking hooks will remain in the same position as in the previous UNLOADING operation (see image C).
- Move the stretcher toward the hooks on the hook’s assembly and press it to the hooks and they will automatically raise the stretcher to the initial position (see image A).
- Using the “EXTRACTION” handle (see image 1), push the upper part of the support into the ambulance.
- Ensure that the stretcher-lift assembly is completely immobilized.

UNLOADING

- Remove the upper part of the support by pressing the “EXTRACTION” handle (see Figure 1) and pulling it out of the ambulance (see Figure A).
- Using the button located on the right side of the upper carriage (see image 2), press (I) to initiate the AUTOMATIC UNLOADING operation; wait until the upper part of the support is fully tilted (see image B).

05 OPERATION


REV. 2026/03

- Immediately afterward, the drive head will begin to descend until it stops (see image C)
- The hooks will automatically release the stretcher; take care to support it during the descent.
- Once the stretcher is on the floor, use the side button (II) to return the upper carriage to its initial resting position (see image D); the hook's assembly will remain in the same position to reduce loading time.
- Using the "EXTRACTION" handle (see image 1), push the upper carriage into the ambulance.
- Ensure that the electric stretcher holder is completely immobilized.

05 OPERATION


REV. 2026/03

LATERAL DISPLACEMENT

Operate the pedal located on the bottom rear part of the support by applying a slight pressure with your foot and guiding the movement by sliding it with your hands.

The total lateral displacement is 300 mm.


05 OPERATION

REV. 2026/03

AUTOMATIC STOP DUE TO OVERPRESSURE DETECTED BY THE FRONT SENSORS

During the loading process, the system stops if overpressure is detected by the front safety sensors (No. 3 in section 05.1)


The procedure to follow after a shutdown triggered by the front sensors is the following:

- Unload the stretcher manually using the manual release bar (11).
- Resolve the situation that caused the shutdown.
- Return the electric stretcher holder PB-440 to the ZERO or REST position by pressing button No. 6 in Figure 05.1, position II.
- Wait for the platform to level out and the carriage to return to its starting position.

STOP VIA EMERGENCY STOP BUTTON (N.7 in figure 05.1)

At any time, it may become necessary to abruptly stop the lift's movement. This can be done by pressing the stop button No. 7 shown in Figure 05.1.

The procedure to follow after a shutdown triggered by the front sensors is the following:

- Unload the stretcher manually using the manual release bar (11).
- Resolve the situation that caused the shutdown.
- Return the electric stretcher holder PB-440 to the ZERO or REST position by pressing button No. 6 in Figure 05.1, position II.
- Wait for the platform to level out and the carriage to return to its starting position.

05 OPERATION

REV. 2026/03

05.5. Electric movements inside the ambulance

Using the internal control inside the ambulance, the following predefined movements can be performed:

ZERO POSITION OR REST POSITION

The zero position, also known as the rest position, is when the stretcher support is completely horizontal and folded in on itself (image below)

TRENDELEMBURG (Protruding / Reclining)

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The movement occurs while the button is pressed (9.I or 9.II). To stop the movement, stop pressing the button.

- Protruding Trendelenburg position:

In this position, the headrest rises diagonally from the position where the patient's head is located; it is activated using button (No. 9-II)

- Trendelenburg decline position:

In this position, the patient's feet are raised diagonally upward; it is activated using button (No. 9-I).

05 OPERATION

REV. 2026/03


05.6. Manual operation of the support

In the event of a power failure, the stretcher support is equipped with a manual operation system.


To operate in manual mode, the stretcher support must be turned off and the carriage must be unloaded. If the stretcher is secured to the carriage, release it using the manual release bar (item 11 on section 5.1).

To operate in manual mode, you must release the actuator from its position. To do so, follow these steps:


1. Remove the retaining clip from the actuator shaft.

2. Release the actuator shaft and let it rest on the lower support.

3. Remove the tap from the back of the motor cover and insert the crank located at the bottom of the bench. Turn the crank until the carriage moves to the top.

05 OPERATION

REV. 2026/03

Once the carriage is at the top, you can operate the stretcher support using the front handles.

Once the entire process is complete, be sure to remove the crank and store it in its designated place (see photo below)


05.7. Electrical diagram

PRINTED CIRCUIT BOARD COMPONENTS


N°	DESCRIPTION OF THE COMPONENTS	N°	DESCRIPTION OF THE COMPONENTS
1	Loading position (Control command)Crank	14	Tilt switch 2
2	Micro-sensor (Loading handle)	15	Tilt switch 1
3	Zero position (Control command)	16	Micro-switch for loading position (common)
4	Micro-sensor (Unloading handle)	17	Micro-switch for loading position
5	Emergency stop	18	Linear Actuator A1 (common)
6	Subsequent signals (common)	19	Linear Actuator A1
7	Mains connection 12V +	20	Motor
8	Mains connection 12V -	21	Motor (common)
9	Trendelemburg declive (Control command)	22	LED (common)
10	Trendelemburg proclive (Control command)	23	LED
11	Upload platform (Control command)	24	Buzzer (siren) (common)
12	Download platform (Control command)	25	Buzzer (siren)
13	Common (Control command)	26	Encoder

05 OPERATION

REV. 2026/03


FEMALE CONNECTOR DIAGRAM

05 OPERATION


REV. 2026/03

MALE CONNECTOR DIAGRAM

CABLE DIAGRAM

-
- A 12V DC power input line (2 wires)
-
- B Control unit connection line (7 wires)
-
- C Actuator A1 (feet) connection line (2 wires)
-
- D Actuator A2 (head) connection line (2 wires)
-
- E Communication line 1 (6 wires)
-
- F Communication line 2 (6 wires)
-
- G Emergency stop line (2 wires)
-
- H Load handle microswitch line (2 wires)
-
- I Unload handle microswitch line (2 wires)
-
- J Load position microswitch line (2 wires)

06 GENERAL MAINTENANCE *REV. 2026/03*

06.1. Cleaning

It is essential to keep the equipment clean to ensure proper use and durability of the set. Thorough cleaning should be carried out periodically, especially in areas exposed to dirt that can accumulate dirt such as corners or holes.

Do not use high pressure cleaning systems or bleach to disinfect the product, they can damage it. Instead, use a hydroalcoholic-based disinfectant, wash it with water and let it dry naturally, do not use direct heat sources to dry.


Non-compliance with cleaning operations can lead to the risk of cross infections due to the presence of secretions and/or residues. During all control and disinfection operations, the operator must wear appropriate personal protective equipment, such as gloves, glasses, etc.

When the cleaning procedure includes a disinfectant product, carefully follow the instructions of the manufacturer of the product used regarding the method of application and contact time. Make sure that you have taken all the necessary precautions to ensure that there are no risks of cross-infection or contamination of patients and operators.

Failure to carry out the correct cleaning routine could increase the risk of cross infection, due to presence of body fluids and/or residuals.
The operator must always wear adequate personal protection such as gloves and mask, etc. during all checking and cleaning procedures.

Make sure to keep these areas free of water and humidity. Take special care not to wet them while washing the set.

06.2. Lubricating

In general, all moving parts must be greased. Our products leave the factory completely greased and lubricated, however, it is possible that over time and use of the product the elements become degreased, either due to loss of lubrication or due to dirt.

Periodically clean and oil affected areas according to manufacturer's specifications. Check for loose, missing or worn parts.

Periodically inspect all moving parts to ensure that components are tight.

06 GENERAL MAINTENANCE *REV. 2026/03*

06.3. Wear areas

Periodic inspection of the different components of the system for signs of wear is a preventative measure that can reduce breakdowns. Check especially the moving parts such as the ramp and the lateral shift pedal, if present.

In the event that a worn component needs to be changed, you should contact our sales department for more information on ordering spare parts and their installation.

The repair or replacement of product components must necessarily be carried out by qualified personnel, that is, who have received training from the Manufacturer, Promeba S.L. and original spare parts should always be used.

Promeba S.L. declines any responsibility for any damage, direct or indirect, resulting from improper use of spare parts and/or any repair intervention carried out by unauthorized persons.

06.4. Maintenance program

Although the frequency of the controls is determined by factors such as legal requirements, the type of use, the frequency of use, the environmental conditions during use and storage, it is placed a recommendation on the maintenance that the equipment must have for correct operation throughout its useful life.

All maintenance and revision activities must be recorded and documented with the corresponding technical intervention reports, and the documentation must be kept for at least 10 years from the end of the product's useful life and must be made available to the competent Authorities and/or from the Manufacturer, when requested.

Maintenance review	Every use	When necessary	Every month	Every year
Disinfect	X		X	X
Clean		X	X	X
Inspect	X	X	X	X
Grease moving parts and the lateral displacement				X
Change the scoop stretcher clips				X
Check the tension of the lateral displacement cables		X		

For both the disinfection and cleaning procedures, point 7.1 of this manual must be followed.

During the inspection procedure, which must be done after each use, the following tasks must be carried out:

- Check that all components are present.
- Check the integrity of the device, that is, that there are no breaks, cracks, holes or cuts.
- Check the wear of the product following section 7.3 of this user guide.
- Check that the moving parts slide correctly.
- Check that the accessories provided are present, intact and that they function correctly.


07 SPARE PARTS

REV. 2026/03

To request a spare part, indicate:

- Number and description of the set in which the necessary part is found (see table)
- Point out the piece in the drawing, or the area where it is located.

Following these instructions, our sales department will provide you with a detailed drawing of the assembly so that you can specify the reference of the required spare part.

N°	COMPONENTS DESCRIPTION
1	TOP PLATFORM
3	ELEVATOR GUIDE SET

N°	COMPONENTS DESCRIPTION
2	BASE PLATFORM

10 LEGAL NOTICES

REV. 2026/03

- This document may contain technical inaccuracies or typographical mistakes.
- Changes are periodically added to the information herein; these changes will be incorporated in new editions of the publication.
- Promeba, S.L. reserves the right to make any modification or improvement in the products described in this publication if it is appropriate.
- Promeba, S.L. may have patents or patent applications which address themes described in this document. The possession of this document does not entitle no license to these patents. The information contained in this document does not affect or change the specifications or product specifications or warranties of Promeba, S.L.
- No part of this document shall operate as express or implied license or compensation under the intellectual property rights of Promeba, S.L. or .third party.
- All the information contained in this document has been obtained in specific environments and it is presented as an illustration. The results obtained in other operating environments may vary.
- Promeba, S.L. can use or distribute the information that the customer provides in whatever way they see fit, without incurring any obligation with the client.

11 PRODUCT WARRANTY *REV. 2026/03*

Promeba, S.L. guarantees that its products have satisfactorily passed all established quality controls, both functional and material. The duration of the guarantee is 2 years from the date of purchase of the product.

This guarantee will be valid only when it is presented with the original invoice or proof of purchase (indicating the date of purchase, model and the name of the distributor) together with the defective product during the period covered by the guarantee. Promeba, S.L. reserves the right not to offer the free warranty service if the indicated documents are not presented or if the information they contain is incomplete or illegible.

1. This warranty will not apply if the model name or serial number of the product has been altered, erased, disappeared or becomes illegible.

2. This guarantee does not cover transportation costs or risks derived from the transportation of your product to and from Promeba, S.L.

3. This warranty does not cover any of the following:

- a) Periodic maintenance and repair or replacement of parts derived from normal wear and tear.
- b) Expendable material (components that are expected to need periodic replacement during the life of the product, such as belts, straps, etc.).
- c) Damages or defects derived from the improper use, operation or treatment of the product and not due to normal use of the product.
- d) Damages derived from:
 - i. Misuse, including:
 - Treatment that results in damages or physical, superficial or appearance changes of the product.
 - Installation, use or storage of the product in a way that does not

respect the instructions described by Promeba, S.L.

- Maintenance of the product in a way that does not respect the instructions of Promeba, S.L. for proper maintenance.

- Installation or use of the product in a way that does not respect the technical or safety regulations of the country where it is used or installed.

ii. Use of components not provided with the product or incorrect installation of accessory parts not previously tested.

iii. States or defects of the system in which the product is used or incorporated with the exception of other products Promeba, S.L. designed for use with the product.

iv. Use of the product with accessories, peripheral units and other products of a type, condition or standards not established by Promeba, S.L.

v. The manufacturer or distributor will be solely responsible for deciding whether to send the parts for repair, or to replace the product in its entirety. In no case will operators be sent to repair or replacement of the product.

Except in the cases mentioned above, Promeba, S.L. will make no warranties in relation to product, performance, accuracy, reliability, or fitness for logic purpose of the equipment or other type. If this exception is not legal or contemplated by current law, Promeba, S.L. will limit or exclude your warranties only to the extent permitted by applicable law.

The only obligation on the part of Promeba, S.L. in connection with this warranty is to repair or replace parts subject to the terms and conditions of this warranty.

Promeba, S.L. is not responsible for loss or damage to products, this warranty or others, including economic loss or non-assessable

11 PRODUCT WARRANTY *REV. 2026/03*

damage; the price paid for the product; loss of profits, income, information, usufruct or use of the product or associated products or indirect, accidental or critical loss or damage.

This clause refers to whether the loss or damage is due to deterioration or inoperability of the associated product due to defects or unavailability of Promeba, S.L., which has caused downtime, loss of user time or business interruption.

In cases where the law prohibits or limits these liability exclusions, Promeba, S.L. will exclude or limit his liability only to the extent permitted by applicable law.

For example, there are countries that prohibit the exclusion or limitation of damages caused by negligence, reckless negligence, willful misconduct, fraud and similar actions. The responsibility of Promeba, S.L. In this guarantee will not exceed, in any case, the price paid for the product, but if the current law allows only limitations of greater responsibilities, these will apply.

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*All rights reserved. Variations can be done whitout notice.
Promeba, S.L. is to be considered not responsible for damages caused by the lack or the wrongness of the information here mentioned.*

SCAN QR CODE

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Ctra C-16 Km 59.5 · 08650 Sallent (Barcelona) · SPAIN

T. 93 837 12 00

promeba@promeba.com · www.promeba.com